English

12. FRACTURE REDUCTION AND ALIGNMENT

The Aligner PHX Humeral Fracture Brace is designed with adjustable humeral uprights, to allow the physician to control bone alignment during the healing process. This alignment should be done by an ORTHOPEDIC SURGEON, with the assistance of x-ray or fluoroscopy, to ensure proper fracture reduction. Adjustment(s) may be necessary to maintain acceptable alignment.

13. FRACTURE REDUCTION PROCESS

The humeral arm has two adjustment pivot points to allow fracture reduction. Either hinge can be used individually or in combination to facilitate proper fracture alignment.

- 1: Top Adjustment
- 2: Bottom Adjustment

The varus/valgus hinge(s) are adjusted using the enclosed Allen wrench. Make sure to tighten the hinges once adjustments have been made.

14. HINGE INSTRUCTIONS

- 1: Flexion Limit Button (ROM)
- 2: Extension Lock Tab
- 3: Extension Limit Button (ROM)
- 4: Locking Degrees of Extension
- 5: Indicator Mark

Depress Extension Button and slide to desired ROM setting. Depress Flexion Button and slide to desired ROM setting.

NOTE: ROM settings automatically lock into place.

16. LOCK HINGE

The Lock Hinge allows the brace to be locked at a specific setting. The hinge can be locked in 10° increments ranging from -10 to 100.

Bend the arm to desired setting and engage lock. This will lock the brace in place independent of ROM settings. Make sure ROM settings do not interfere with desired lock setting.

THERE ARE THREE WAYS TO LOCK THE ELBOW HINGES:

- Use extension lock tab.
- ROM limits. Set ROM hinge to desired settings by pushing in the gray tabs and sliding them along the hinge to desired position.
- Locking the ROM limits.

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INDICATIONS: For fractures of the mid-shaft or distal third humerus. Recommend placement 5 to 7 days post injury or under discretion of the treating physician.

CONTRAINDICATIONS: Intra-articular distal humeral fractures. Not appropriate for treatment of open fractures.

CAUTION: This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is not intended to prevent

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
injuries or reduce or eliminate the risk of reinjury. **WARNING:** Remove or adjust only under direction of your physician.

WARNING: If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.

English

INDICATIONS: For non-operative distal to mid-shaft humeral fractures.

CLINICIAN FITTING INSTRUCTIONS:

1. APPLY THE SHOULDER HARNESS
Apply to the unaffected shoulder and tighten to fit snugly under the arm with the back pad centered along the spine. Orient the pad as shown.

2. FITTING THE ALIGNER PHX BRACE
Hold the brace next to the patient's arm, and align the hinge with the elbow axis. The deltoïd cuff should cover the shoulder joint. The forearm should extend from the elbow axis to the wrist.

3. ADJUSTING THE HUMERAL AND FOREARM SECTION

(A) To adjust, press the silver button and slide the humeral section to the desired length. Release the button to click into place. The deltoïd cuff should be on top of the shoulder.

(B) Press the silver button on the forearm section of the brace and slide to the desired length. Release the button to click into place.

4. LOCK THE ELBOW HINGE
Set the hinge at 90 degrees and then lock.

5. LOOSEN ALL THE STRAPS
Loosen straps on the brace so the arm can be placed between the medial and lateral humeral shells.

6. PLACE THE BRACE ON THE ARM
The forearm should rest in the forearm cuffs and the elbow should be aligned with the center of the hinge.

7. FITTING THE HUMERAL SECTION
Begin fitting by positioning the lateral humeral cuff against the outer arm and position the inner shell against the medial side of the arm

so that the epicondylar flare lines up with the medial epicondyle. This will create the clamshell around the humerus. If too long, the cuff can be trimmed; however, leave enough length so the shells overlap slightly. The upright should be centered on the arm.


8. FITTING THE HUMERAL CUFFS
While holding the cuffs in place, wrap the strap ends around the arm and through the D-rings. Make sure no skin is pinched between the cuffs as you pull the straps to the desired tension. Secure the hook and loop to the straps adjusting both anterior straps (1) first, then both posterior straps (2).

9. FITTING THE FOREARM SECTION
Mold the forearm cuffs to the patient by squeezing the moldable aluminum shell. The cuff should be snug but not uncomfortable. Mold both of the cuffs to the patient's forearm.

10. SECURE STRAPS ON THE FOREARM
To secure the straps, first engage the quick release buckles. Disengage the hook and loop on the straps and pull straps to the desired tension. Secure the strap with the hook and loop attachments. Trim off any excess strap material as needed.

11. ENGAGE THE SHOULDER HARNESS
The harness straps form a figure eight around the shoulders. With the back pad positioned in the middle of the back, feed straps through the D-rings on the humeral cuff and then back through the buckle on each strap and adjust the tension, making sure the back pad stays in the center of the back.

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