



17AF10=*

DE Gebrauchsanweisung	3
EN Instructions for use	4
FR Instructions d'utilisation.....	6
IT Istruzioni per l'uso	7
ES Instrucciones de uso	9
PT Manual de utilização.....	10
NL Gebruiksaanwijzing.....	12
SV Bruksanvisning	13
DA Brugsanvisning	15
NO Bruksanvisning	16
FI Käyttöohje.....	18
PL Instrukcja użytkowania.....	19
HU Használati utasítás	21

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INFORMATION

Last update: 2017-07-03

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

These instructions for use provide you with important information on the processing of the 17AF10=* ankle joints.

2 Intended use

2.1 Indications for use

The product is intended **exclusively** for use in pairs, in orthotic fittings of the lower limbs and **exclusively** for contact with intact skin. The orthosis must be used in accordance with the indications.

2.2 Indications

- Dorsiflexor weakness

Indications must be determined by the physician.

3 Safety

3.1 Explanation of warning symbols

⚠ CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Incorrect alignment, assembly or adjustment

Injury and damage to product components

- ▶ Observe the alignment, assembly and adjustment instructions.
- ▶ Work on the product may only be carried out by qualified personnel.

CAUTION

Use of the product in case of known allergies to chrome or nickel

Allergic reaction

- ▶ Avoid skin contact with the product in case of known allergies to chrome or nickel.

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

4 Preparation for use

4.1 Environmental conditions

Permitted: Temperature range for use -10°C to $+60^{\circ}\text{C}$; relative humidity 0 % to 90 %, non-condensing.

Not permitted: Mechanical vibrations, impacts, perspiration, urine, fresh water, salt water, acids, dust, sand, highly hygroscopic particles (such as talcum).

4.2 Processing

For further processing instructions, see the 646T523=* technical information.

- 1) Place the dummy on the joint pivot point of the plaster model and secure it with the nails (see fig. 1).
- 2) Vacuum-form or laminate the orthosis.
- 3) Remove the orthosis and dummies from the plaster model after it has hardened/cooled.
- 4) Separate the foot and calf component with a saw and deburr the edges (see fig. 2).
- 5) For pre-flexed joints make a U-shaped cutout, for straight joints make a V-shaped cutout at the front sides of the joint adapters (see fig. 3).
- 6) Drill the mounting holes for the joints (diameter 4.5 mm).
- 7) Insert the joints and attach them to the orthosis with the screws (see fig. 4) (torque value: 0.5 Nm).
- 8) Cover both joints with the padding.

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the

manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH