sherpa3D

powered by OssKin

Digitally bespoke 3D-printed
multi-material foot orthotics

Corrective Elements Guide

Table of Contents

1. Elements & Sinks

1. Metatarsal Pad
2. Metatarsal Bar
3. Met. Head Cushion (1-5)
4. Met. Head Cut Out (any of 1-5)
5. 1st Digit Offload Pad
6. Kinetic Wedge
7. Toe Crest Pad (1-5)
8. Toe Crest Pad (2-5)
9. Toe Crest Pad (1)
10. Cluffy Wedge
11. Heel sink

2. Forefoot Options

1. Internal Wedge
2. External Wedge

3. Midfoot Options

1. Additional Medial Arch Height
2. Anterior Arch Extension
3. Lateral flange
4. Arch rigidity

4. Rearfoot Options

1. Heel cup depth
2. Heel lift
3. Internal post / wedge
4. External post / wedge

1. Elements & Sinks

The following topographical elements are additions or modifications to the plantar surface of the orthotic shell. Each element targets specific anatomical structures to redistribute plantar pressure, accommodate deformity, or promote optimal gait mechanics.

Overview - All Elements & Sinks

Detailed - Description for each Element & Sink

Metatarsal Pad

Function:

Promotes restoration of the transverse metatarsal arch by elevating the metatarsal shafts proximal to the heads. This redistributes plantar pressure away from the central metatarsal heads, effectively reducing compression on interdigital nerves.

Metatarsal Bar

Function:

Transfers loading forces away from the entire metatarsal head row (M1–M5) by providing a supportive ridge immediately proximal to the metatarsophalangeal joints (MPJs). Effective for diffuse metatarsalgia and reducing peak pressure during the terminal stance phase of gait.

1. Elements & Sinks

Met. Head Cushion (1-5)

Function:

Provides localized shock absorption and pressure attenuation beneath metatarsal heads via low-durometer padding. It accommodates fixed plantarflexed metatarsals or localized fat pad atrophy without changing the underlying skeletal alignment.

Met. Head Cut Out (any of 1-5)

Function:

Creates a localized drop (sink) within the Met. Head Cushion directly beneath the targeted metatarsal head. This maximizes offloading for prominent, rigid, or ulcerated metatarsal heads by redirecting ground reaction forces to adjacent structures.

1st Digit Offload Pad

Function:

Topological elevation of all digits except the 1st digit from the metatarsal heads to the toes. Reduces localized pressure during the propulsion phase; frequently utilized for diabetic ulcer prevention or hallux rigidus management.

Kinetic Wedge

Function:

Features a localized structural drop-out zone specifically under the first metatarsal head, combined with extension support under digits 1-5. This permits the first metatarsal to plantarflex naturally, facilitating sagittal plane progression.

Toe Crest Pads (1-5; 2-5; 1)

Function:

Fills the sub-digital sulcus space to reduce apical pressure on clawed, hammer, or mallet toes. By supporting the proximal phalanges, it relieves tensile stress on the long flexor tendons. Available in three configurations: all toes (1-5), lesser toes only (2-5), or hallux only (1).

1. Elements & Sinks

Cluffy Wedge

Function:

A localized small wedge placed directly under the hallux proximal phalanx, promoting pre-extension of the first MPJ. This pre-loading assists in overcoming functional hallux limitus, encouraging timelier activation of the plantar fascia windlass mechanism during late stance.

Heel Sink

Function:

A localized pocket filled with a low-durometer insert positioned within the center of the heel cup. Attenuates peak ground reaction forces at initial contact (heel strike), accommodating calcaneal spurs or localized plantar fasciitis.

2. Forefoot Options

Forefoot wedging modifies the frontal-plane orientation of the forefoot relative to the rearfoot, correcting deformity or optimizing loading across the metatarsal heads. Available in up to +3 mm increments.

Overview - Forefoot Options


Figure 1 — Forefoot Options: Internal Wedge (supination) and External Wedge (pronation), front view with +3 mm reference.

Detailed - Description for each Forefoot Option

Internal (Supinating) Wedge

Function:

Controls premature visual midfoot collapse by increasing supination, or supports a forefoot varus deformity. The medial forefoot edge is elevated, shifting the center of pressure laterally and reducing excessive pronation forces during the loading response.

External (Pronating) Wedge

Function:

Optimizes lateral column loading by increasing pronation, or supports a forefoot valgus deformity. The lateral forefoot edge is elevated, shifting the center of pressure medially and accommodating a rigid plantarflexed first ray.

3. Midfoot Options


Midfoot modifications alter the arch geometry and rigidity of the orthotic shell to control pronation, distribute plantar forces, or provide structural containment.

Detailed - Description for each Midfoot Option

Additional Medial Arch Height

Function:

Increases the vertical profile of the shell to maximize direct contact with the medial longitudinal arch. By filling the void under the midfoot, it expands the weight-bearing surface area, reducing peak pressures in the forefoot and rearfoot. This directly resists excessive foot pronation, stabilizing the medial column against structural collapse..

Anterior Arch Extension

Function:


Extends the longitudinal arch apex further forward toward the metatarsal heads. This increases total surface contact area under the midfoot, distributing forces along the entire plantar fascia complex.

Lateral Flange

Function:

A rigid vertical extension along the lateral border of the shell. Acts as a structural physical block to counteract severe lateral translation, subluxation, or excessive supination forces.

Arch Rigidity


Function:

Modifies the physical flex-modulus and resistance to deformation of the midfoot shell structure. Increasing rigidity controls high-mass loading or rigid pronation; lowering rigidity permits functional midtarsal mobility and dynamic shock absorption.

4. Rearfoot Options

Rearfoot modifications influence calcaneal containment, limb-length accommodation, and subtalar joint axis alignment — providing proximal control of the entire kinetic chain.


Detailed - Description for each Midfoot Option

Heel Cup Depth

Function:


Controls the vertical containment wall around the calcaneum. Available in three profiles: 8 mm (standard), 13 mm (moderate), and 18 mm (deep). Higher walls maximize calcaneal control, reducing calcaneal eversion/inversion deviations relative to the lower leg.

Heel Lift

Function:


Elevates the rearfoot platform to accommodate limb length discrepancies (LLD) or to reduce tension on the Achilles tendon and posterior chain musculature. Available up to +5 mm (default: 0 mm).

Internal Post / Wedge

Function:

A structural stabilizing addition designed into the plantar surface of the shell geometry (intrinsic varus post). Modifies the subtalar joint axis alignment by shifting the center of pressure medially, directly at the foot-orthotic interface.

External Post / Wedge

Function:

A structural stabilizing addition designed into the plantar surface of the shell geometry (intrinsic valgus post). Modifies the subtalar joint axis alignment by shifting the center of pressure laterally, directly at the foot-orthotic interface.

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2. Options avant-pied *(Forefoot Options)*

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