

SELECTION[®] ULNAR DEVIATION SPLINT

Improve finger alignment to enhance hand function and functional grip.



- Adjustable fastening to optimize the directional pull
- Low profile, breathable fabric
- Shaped to the contours of the hand
- Trimmable with scissors
- Machine washable

SELECTION®

ULNAR DEVIATION



SELECTION® Ulnar Deviation Splint is intended to improve finger alignment and enhance hand function and functional grip.

This hand splint is exclusively designed using our proprietary three layer Selection® fabric. With its low profile design, it meets the contours of the hand, without impeding function. The fabric can be trimmed without fraying, and therefore allows further customization if required. The easy fastening surface allows the hook to be placed anywhere on the spelling of optimize the directional pull needed. The orthosis is designed so that the material extends far to the ulnar side at the 5th MCP to support the best possible direction of pull and distribute the pressure.

The Ulnar Deviation splint supports the arches of the hand while providing realignment and stabilization, which can contribute to the management of pain.



Intended use

For patients with inflammatory and degenerative diseases, such as rheumatoid arthritis, that has resulted in ulnar deviation. the Selection® Ulnar Deviation splint can also be used for hypermobile joints and weak muscles to provide stability. Not intended for use with fixed contractures in the hand.



Item No. Left	Item No. Right	Size	Circumference MCP 2-5
35218 1010	35218 2010	X-Small	14 - 17 cm
35218 1011	35218 2011	Small	17 - 20 cm
35218 1012	35218 2012	Medium	20 - 23 cm
35218 1013	35218 2013	Large	23 - 27 cm