



# WALK A MILE IN THEIR SHOES. (OR 2000)

Having undergone over 2 million cycles of testing, the M2 AFO™ has proven to be durable for extended use to facilitate a normal gait. The lightweight and comfortable design offers ideal support with optimal flex for toeoff.








**Carbon Fiber**  
Creates dynamic forces for comfortable walking.



**Trimmable Footplate**  
Can be customized to fit most patient anatomies and shoes.



**Soft Padding**  
Provides comfort for all-day wear.



**Medium Flex**  
Fits most patient needs for enhanced performance.

POSTERIOR M2 AFO - Code L1951 Suggested			ANTERIOR M2 AFO - Code L1932 Suggested	
SIZE	LENGTH	HEIGHT	LENGTH	HEIGHT
S	10.04 in (255 mm)	13.78 in (350 mm)	10.04 in (255 mm)	15.35 in (390 mm)
M	10.63 in (270 mm)	14.17 in (360 mm)	10.63 in (270 mm)	15.75 in (400 mm)
L	11.42 in (290 mm)	14.57 in (370 mm)	11.42 in (290 mm)	16.14 in (410 mm)
XL	11.61 in (295 mm)	14.96 in (380 mm)		

POSTERIOR	RIGHT		LEFT	
	Small - Code L1951 Suggested	640112	Small - Code L1951 Suggested	640212
	Medium - Code L1951 Suggested	640113	Medium - Code L1951 Suggested	640213
	Large - Code L1951 Suggested	640114	Large - Code L1951 Suggested	640214
	XLarge - Code L1951 Suggested	640115	XLarge - Code L1951 Suggested	640215

ANTERIOR	Small - Code L1932 Suggested	640102	Small - Code L1932 Suggested	640202
	Medium - Code L1932 Suggested	640103	Medium - Code L1932 Suggested	640203
	Large - Code L1932 Suggested	640104	Large - Code L1932 Suggested	640204

ACCESSORIES	Anterior Replacement Pads, Right	640101
	Anterior Replacement Pads, Left	640201
	Posterior Replacement Pads, Universal	640111