



Buttress Support Soft Knee Brace

Sizes: xs, s, m, l, xl, xxl

The Buttress Support's 1/4" foam buttress can be positioned inside the neoprene pocket to offer stabilization from virtually any position.



Common Examples of Use

- Patella tendonitis
- Patella subluxation
- Chronic patellar maltracking
- Mild chondromalacia

Product Features

- Straps above and below knee
- Medial and lateral stays help prevent migration
- Also available with an open back

Category: Patellofemoral.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
0721X	Buttress Support, Open Back

Product Sizing Chart

PTO / Patella Supports (All)			
Size #	Size	Thigh Circumference	
X = 1	XS	12" – 15"	(30 – 38 cm)
X = 2	S	15" – 18"	(38 – 46 cm)
X = 3	M	18" – 21"	(46 – 53 cm)
X = 4	L	21" – 24"	(53 – 61 cm)
X = 5	XL	24" – 27"	(61 – 69 cm)
X = 6	XXL	27" – 30"	(69 – 79 cm)
<i>Circumference taken 6" (15 cm) above mid patella</i>			