


ANKLE for SINGLE AXIS FOOT

FOR ENDOSKELETAL PROSTHESIS

• Ankle for 1B02 Single Axis Foot

Item N°.	Max patient's Weight	Material	Weight	Upper connection
1D12	100 kg	Alu	250 g	Tube Ø 30 mm



• Selection of False Malleoli

Item N°.	Description	Material	Weight
1E1301	For foot from 22 to 24 cm	PU	25 g
1E1401	For foot from 25 to 27 cm	PU	30 g
1E1501	For foot from 28 to 29 cm	PU	35 g



1E1301

FOR EXOSKELETAL PROSTHESIS

• Malleolar Joint and Leg Base

	Item N°.	Description	Max patient's Weight	Material	Weight
Joint	1D02	For foot from 22 to 24 cm	100 kg	Alu	250 g
	1D03	For foot from 25 to 27 cm	100 kg	Alu	260 g
	1D04	For foot from 28 to 29 cm	100 kg	Alu	270 g
Leg base	1D22	For foot from 22 to 24 cm	100 kg	Wood	580 g
	1D23	For foot from 25 to 27 cm	100 kg	Wood	590 g
	1D24	For foot from 28 to 29 cm	100 kg	Wood	600 g



1D02



1D23

RELATED PRODUCTS

• Selection of Rear Bumpers for 1B02 Foot



Item N°.	Hardness	Thickness	Color	Material	Weight
EE013	X-soft	50 mm	Brown	Elastomer	17 g
EE007	Soft	50 mm	Beige	Elastomer	20 g
EE008	Medium	50 mm	Black	Elastomer	20 g
EE009	Hard	50 mm	Gey	Elastomer	20 g

• Selection of Front Bearings



Item N°.	Thickness	Description	Material	Weight
EE011	5 mm	For 10 mm heel	Elastomer	2 g
EE001	6 mm	For 14 mm heel	Elastomer	3 g
EE002	7 mm	For 19 mm heel	Elastomer	4 g
EE003	9 mm	For 24 mm heel	Elastomer	5 g

• 5F055 Key for Malleolar Nut