

CONNECT[®] TF

Expanded Access

Connect TF helps address the challenges that lead to fewer than 50% of low active transfemoral amputees using prostheses.

The prosthetist can easily make adjustments in response to limb changes. Patients can comfortably wear Connect TF, as well as easily don and doff, while seated.

- ① Prosthetist can easily adjust height, volume and shape.
- ② Soft proximal brim and conformable shells for patient comfort, whether sitting, standing, or walking.
- ③ On/off tensioning handle expands circumference approximately 1" (2.5cm) for easy donning/doffing.

PRODUCT HIGHLIGHTS

- Prosthetist can quickly adjust Connect TF to adapt to limb changes.
- Easy to put on and to use by low active or new amputees.
- Custom fit and delivery in clinic, hospital, rehab, home or nursing facility.

Exclusively available in the Össur Legs program.



CONNECT® TF

POTENTIAL PATIENTS

- Amputee who can't don traditional socket while standing
- New amputee with rapidly changing limb volume
- Existing amputee that rejects traditional socket over comfort or donning

Patient Information	
Amputation Level	Transfemoral
Impact Level	Low
Patient Weight Limit	125kg (275lbs)
Socket Information	
Sizes	Medium-Short, Medium-Standard, Large-Short, Large-Standard
Weight of Socket	1.5 kg (3.3lbs)
Liner Selection	Iceross Transfemoral Locking
Limited Warranty	24 Months

ÖSSUR LEGS

When you purchase the Connect TF socket in an Össur Leg package, you can select from a wide choice of clinically appropriate options below.

Description	
Suspension	Locking
Socket	Connect TF with clutch lock and 4-hole pattern
Liners	Iceross Transfemoral Locking
Knee	Balance Knee OFM1 or Balance Knee OFM2
Adapters	Clinically needed Össur adapters for the amputee (including Titanium)
Foot & foot cover	Balance Foot S, K2 Sensation, Balance Foot J or no foot
Options	Balance Knee OFM1 with IKF adapter Total Knee 1900 Locking Knee D/P Flexion hydraulic ankle Balance Foot S Torsion

Have you tried Össur Legs?

Ask your Össur Area Manager for an example of how it benefits you.



FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM
CANADA (800) 663-5982
WWW.OSSUR.CA

