

MEASUREMENT CHART

DDS 500

Lumbar Spinal-Air
Decompression LSO Brace



DDS Double

Lumbar Spinal-Air
Decompression LSO Brace



DDS 500 & DDS Double Sizing Chart			
DDS 500 Brace Size	DDS Double Brace Size	Torso Measurement	Panel Size
S	S	26" -28"	#1/ Small
M	M	29" -32"	#2/ Medium
L	L	33" -35"	
XL	XL	36" -38"	
2XL	2XL	39" -41"	#3/ Large
3XL	3XL	42" -44"	
4XL	4XL	45" -48"	#4/ X-Large
5XL	5XL	49" -51"	
6XL	n/a	52" -55"	#5/ 2X-Large
7XL	n/a	56" -59"	

The measurement needs to be taken at the most protruding part of the patients "torso", which is usually level with the patients naval.



DDS 8" Extension

Cervitrac™

Cervical Traction Collar



DDS CERVITRAC Sizing Chart



Size	Circumference
SMALL	11.5" - 14.5"
LARGE	14.5" - 17.5"

OA PRO™

Decompression Knee Brace



Size	Thigh (6" above the center of the knee)	Center of Knee	Calf (6" below the center of the knee)
S	15.5"~18.5"	12"~14"	12"~14"
M	18.5"~21"	14"~15"	14"~16"
L	21"~23.5"	15"~17"	16"~18"
XL	23.5"~26.5"	17"~19"	18"~20"



As shown in the figure above, measure the circumference of both the thigh and the calf 6" away from the center of the knee.