



EVENUP SHOE LEVEL

- The Evenup™ offers two thicknesses: 1/2 inch and 1 inch
- Easily attaches to the shoe of the non-affected limb for improved balance, gait and reduced body strain
- Specially designed rubber, hook and loop straps hold the Evenup in place
- Includes a 2 layer non-skid rubber outsole
- Eliminates lower back pain and hip pain often associated with use of a walking boot
- Fits right or left foot
- In the event your shoe sizing falls between two sizes, it is recommended to size down

Small (T15022) fits shoe sole measurement between 10" and 11-1/4" and shoe sizes Women's 5 to 8-1/2, Men's 6 to 8

Medium (T15024) fits shoe sole measurement between 11-1/2" and 12-3/4" and shoe sizes Women's 9 to 11, Men's 8.5 to 10.5

Large (T15026) fits shoe sole measurement 13" to 14-1/2" and shoe sizes Women's 11.5 to 13, Men's 11 to 13

FRACTURE SOCK

- Designed to reduce skin irritation while wearing an ankle walker
- 100% synthetic knit allows for superior moisture absorption
- Improved cross stretch provides better compression; helps decrease the risk of friction and lessens the risk of the sock bunching
- The overall superior knit quality of the sock creates a far more comfortable barrier for the patient's skin
- Designed for all ankle walkers
- Sold individually
- 2 lengths available: high or low top
- Suggested HCPCS Code: L2840



<i>Length</i>	<i>Part #</i>
High Top	T15000
Low Top	T15001