

Alignment Components

Universal Coupling Ring

FILLAUER


- When disconnect and rotation are required

124105 Universal Coupling Ring, with Top Plate


Includes:

- A. 885522** Washer, 2.125"x.070",
- B. 124103** Universal Coupling Ring
- C. 882612** Set Screw Cone Point, M6 x 1 x 12mm (3 Each)

Caution: When vacuum forming this component, use the drape forming method with a seam, not the blister forming method (frame and platen)



Bottom View



Universal Coupling Ring with Screws


FILLAUER

- When rotation and disconnect is required


124106 Adapter Ring, with Screws, Anodized

Includes:

- A. 124106A** Adapter Ring, Anodized
- B. 880288** Socket Cap Screw, M6 x 1 x 30mm Flat Head (4 Each)
- C. 882612** Set Screw Cone Point, M6 x 1 x 12mm (3 Each)



Bottom View



Alignment Components

A-P/M-L Slide Unit with Rotation and Pyramid

FILLAUER


124140 A-P and M-L Slide Unit, with Top Plate Assembly and Pyramid

Includes:

- A. 882511** Socket Head Cap Screw, *M5 x .8 x 10mm*
- B. 124155** Top Plate Assembly
- C. 882416** Socket Head Set Screw, *M4 x .7 x 12mm*
- D. 124133** Gib, *Hardened Steel*
- E. 882408** Socket Head Cap Screw, *M4 x .7 x 8mm*
- F. 124110** Center Plate
- G. 124150** Dome Plate Assembly



Bottom View



A-P/M-L Slide Unit with Pyramid and Socket Attachment Plate

FILLAUER


124141 A-P and M-L Slide Unit, with Pyramid and Socket Attachment Plate

Includes:

- A. 880136** Flat Head Undercut Screw, *10-32 x 5/16"*
- B. 124164** Socket Attachment Plate
- C. 882416** Socket Head Set Screw, *M4 x .7 x 12mm*
- D. 124133** Gib, *Hardened Steel*
- E. 882408** Socket Head Cap Screw, *M4 x .7 x 8mm*
- F. 124110** Center Plate
- G. 124150** Dome Plate Assembly



Bottom View



Alignment Components


A-P/M-L Slide Unit with Rotation and Socket Attachment Plate

FILLAUER

124142 A-P and M-L Slide Unit, with Top Plate Assembly and Socket Attachment

Includes:

- A. 882511** Socket Head Cap Screw, $M5 \times .8 \times 10mm$
- B. 124155** Top Plate Assembly
- C. 882416** Socket Head Set Screw, $M4 \times .7 \times 12mm$
- D. 124133** Gib, *Hardened Steel*
- E. 124110** Center Plate
- F. 124164** Socket Attachment Plate
- G. 880136** Flat Head Undercut Screw, $10-32 \times 5/16"$




A-P/M-L Slide Unit with Pyramid Receiver and Rotation

FILLAUER

124143 A-P and M-L Slide Unit, with Pyramid Receiver, Top Plate Assembly

Includes:

- A. 882511** Socket Head Cap Screw, $M5 \times .8 \times 10mm$
- B. 124155** Top Plate Assembly
- C. 882416** Socket Head Set Screw, $M4 \times .7 \times 12mm$
- D. 124133** Gib, *Hardened Steel*
- E. 882408** Socket Head Cap Screw, $M4 \times .7 \times 8mm$
- F. 124110** Center Plate
- G. 124190** Pyramid Receiver Plate
- H. 880290** Socket Set Flat Point, $M8 \times 1.25 \times 16mm$



Alignment Components


A-P/M-L Slide Unit with Two Socket Attachment Plates

FILLAUER


124144 A-P and M-L Slide Unit, *with Two Socket Attachment Plates*

Includes:

- A. 880136** Flat Head Undercut Screw, 10-32 x 5/16"
- B. 124164** Socket Attachment Plate
- C. 882416** Socket Head Set Screw, M4 x .7 x 12mm
- D. 124133** Gib, *Hardened Steel*
- E. 882408** Socket Head Cap Screw, M4 x .7 x 8mm
- F. 124110** Center Plate



Bottom View




Rotatable Socket Adapter with Pyramid Receiver

FILLAUER

124114 Rotatable Socket Adapter Assembly, *with Pyramid Receiver*

Includes:

- A. 124113** Rotatable Adapter
- B. 880290** Socket Head Set Screw, M8 x 1.25 x 16mm



LOWER EXTREMITY PROSTHETICS

Alignment Screws

EAS Metric Screws

FILLAUER

SCREWS	PRODUCT NO.	SIZE	HEAD TYPE	WRENCH SIZE	END USE
	882408 882511 882513 882520 882613 882620 880289 882418	M4x.7x8 M5x.8x10 M5x.8x16 M5x.8x20 M6x1x12 M6x1x20 M8x1.25x50 M4x.7x14	Socket Cap Socket Cap Socket Cap Socket Cap Socket Cap Socket Cap Socket Cap Socket Cap	3mm 4mm 4mm 4mm 5mm 5mm 6mm 3mm	Slide Stop–Dome Plate Slide Stop–Top Plate Pylon Tube Clamp Aluminum Pylon Tube Clamp Titanium Pyramid Receiver Upper Pylon Tube Clamp SACH Foot Bolt Clamp Bolt
	880280 880282 880285 880288 882524	M6x1x16 M6x1x20 M6x1x25 M6x1x30 M5x.8x20	Socket Cap/Flat Socket Cap/Flat Socket Cap/Flat Socket Cap/Flat Socket Cap/Flat	4mm 4mm 4mm 4mm 3mm	Socket Attachment Plate Socket Attachment Plate Socket Attachment Plate Socket Attachment Plate Socket Attachment Plate
	882506	M5x.8x8	Socket Button	3mm	Slide Stop–Counter Sunk Plate
	882612	M6x1x12	Socket Set/Cone Point	3mm	Adapter Ring and Universal Coupling Rings
	882416	M4x.7x12	Socket Set/Cup Point	3mm	Slide Clamps
	880290 882608	M8x1.25x16 M6x1x12	Socket Set/Flat Point Socket Set/Flat Point	4mm 3mm	Pyramid Adjustment Pyramid Adjustment

Alignment Screws

HOSMER

- Used with Spectrum Alignment System

60248	Set Screw, M8 x 12mm
60598	Set Screw, M5 x 6mm, with Nylock
60329	Flat Head Screw, M6 x 20mm
60610	Pan Head Screw, Phillips, 6-32 x 25 Large
56039	Socket Head Screw, M6 x 25mm
60862	Short Screw Slide Shift
60957	Screw Slide Shift

Metric Wrenches

124165M	Metric Wrench Set,
	Hex Size Range 1.5 - 5.0mm
	Tool Type Hex End Keys
	Fastener Type Standard Hex
	Measurement Type Metric
	Handle Type L-Key Long Arm
	Container Type Plastic Holder
	Number of Pieces 6
	Hex Size (mm) 1.5, 2, 2.5, 3, 4, 5
	Material Protanium