

Hero Arm Order Form


Prosthetist Name:		Clinic Number and Name:	
Contact Email:		Opportunity Number or Patient ID	
Purchase Order Number:			

Is your patient under 18? Yes No


Chirality

Left Right


Measurement Affected Arm



Measurement Sound Arm



Measurement Palm Width



Internal or Battery position?

*To be eligible for an Internal battery arm, there must be at least 20.5cm between the end of the residual and the thumb tip

- Internal Battery
- External Battery

Hand Size?

*The small wrist is unavailable with an Internal battery configuration.

- Small Hand, Small Wrist
- Small Hand, Large Wrist
- Medium Hand
- Large Hand

Additional Options

Add a pull through hole?

*Please note, depending on the patient's residual limb, we may need to add additional length to accommodate a pull-through. We will notify you during the design process if this is the case.

Configure for Single Site?

Medial

Lateral

*This can be changed later via the sidekick app. If yes, please indicate which muscle site should be active.

Add 1.5 cm foreshortening?

*Would you like us to include 1.5 cm foreshortening? As the Hero Arm can look longer than the patient's sound arm.

Add sensitivity adjuster tool?

*Please select if you would like to include a free tool to adjust the sensitivity of the EMG sensors with this order - instructions provided. This tool is for clinical purposes only and must **not** be used by patients.

Bi-Lateral Button?

*Replaces the standard release button with a larger disk shaped button, designed to make moving the hand position easier for bi-lateral users. **This will replace the standard button.** The wrist will need to be rotated to attach or remove arm covers.

Hero Arm Cover Choice*?

*Covers may come at an additional cost, please check your pricelist.

Please select **one** option for the texture for the covers:

Plain

Vortex

Harlequin

AzDeco

Aero

Sentinel

Please select **one** option for the colour for the covers:

Black

White

Silver

Vivid Pink

Orange

Blue

Yellow

Green

Red

Blue Grey

Lilac

Pale Pink

Alternatively the below Disney covers are available for **under 18 only**, they must be ordered with the arm (not available to order at a later date).

Iron Man

Frozen

BB-8

R2-D2

Black Panther

Notes:

Hero Arm Cast/Scan Information

Cast Type

Please let us know if you'll be sending a scan file, or a physical cast. Please note, we cannot work from a check socket or direct scan of a patient's limb.

- Digital Scan Physical Cast

Cast/Scan Notes - Please include any important details regarding your cast/scan here:

Please contact your sales rep for the address to ship your cast to.

PLEASE NOTE: Casts will be destroyed after scanning. If you require the cast to be returned, please indicate this below. You will need to arrange a collection once your cast has been scanned.

Cast handling preference:

- Destroy cast after scanning
 Return cast (I will arrange and cover collection costs)

Please submit this completed form along with a copy of your purchase order to support@openbionics.com.