



Hero RGD/PRO Order Form

Prosthetist Name:		Clinic Number and Name:	
Contact Email:		Opportunity Number or Patient ID	
Purchase Order Number:			

Hand Selection:

- Hero PRO Hero RGD

Chirality

- Left
- Right

Size

- Medium Large

- Add a Socket Component Kit?

Please only select this option if you are not also ordering a Hero Flex and would like to make your own socket.

Includes : 1 x Falcon wrist, 2 x Falcon lamination ring, 4 x MyoRays (for inserting MyoPods into socket) and lamination cap.

Please submit this completed form along with a copy of your purchase order to support@openbionics.com.



Hero Flex Order Form

Prosthetist Name:		Clinic Number and Name:	
Contact Email:		Opportunity Number or Patient ID	
Purchase Order Number:			

Specification

Wrist Type



Falcon Wrist



USMC Wrist



Omega Wrist*

*The **Omega Wrist** comes at an additional cost, please refer to your price list or speak to your sales representative for more information.

*The Omega Wrist is **not** directly compatible with Hero RGD/Hero PRO or Activity Attachments, these products are supplied with USMC style inserts.


Chirality

Left Right

Affected Arm

For guidance on any of these measurements, please refer to our Clinician Guidelines, or contact our configuration team at configuration@openbionics.com.

A - Olecranon to Residuum End (cm)



Sound Arm

E - Olecranon to Terminal Device (cm)

This is the distance between the elbow bone and the base of the terminal device (end of Hero-Flex).



There must be a minimum difference in length between measurements A and E of 2.8mm, or select 'Make as short as possible'.

OR

Make as short as possible

Additional Options

Pull-through hole?

Yes No

Please note, depending on the patient's residual limb, we may need to add additional length to accommodate a pull-through. We will notify you during the design process if this is the case.

Add Kidney Condyles?

Yes No

Kidney condyles attach to the inside of the frame between the frame and the socket, to aid suspension. *This is not available for the short version of the Hero Flex.

MyoBand?

Yes No

Add a MyoBand to the order - this allows you to place MyoPods outside of the socket for forearm suspended Flex.

Notes - Please include any important details regarding the arm configuration here:

Cover Selection*

*Covers are optional and may come at an additional cost, please check your pricelist.

Texture

Plain Vortex Harlequin AzDeco

AND

Colour

Black White Silver Red Pale Pink
 Green Blue Orange Lilac Blue Grey Vivid Pink

OR

Special Covers*

*These come at an additional cost

Mr Beast Call of Duty

Hero Flex Cast/Scan Information

Cast Type

Please let us know if you'll be sending a scan file, or a physical cast. Please note, we cannot work from a check socket or direct scan of a patient's limb.

- Digital Scan Physical Cast

Cast/Scan Notes - Please include any important details regarding your cast/scan here:

Please contact your sales rep for the address to ship your cast to.

PLEASE NOTE: Casts will be destroyed after scanning. If you require the cast to be returned, please indicate this below. You will need to arrange a collection once your cast has been scanned.

Cast handling preference:

- Destroy cast after scanning
 Return cast (I will arrange and cover collection costs)