

Icelock 850 Hybrid

& Iceross Seal-In X Locking TF

Built for comfort and stability, the Icelock[®] 850 Hybrid Lock and Seal-In[®] X Locking TF liner brings the benefits of advanced vacuum suspension of the Seal-In system to users with transfemoral amputations of all activity levels who prefer or rely on locking or lanyard suspensions.

The system combines the benefits of elevated vacuum, such as reduced pistoning, volume fluctuations, improved balance, and reduced risk of falls, contributing to comfort in the socket.

The Quick Connect design allows for donning while seating and simple release with the push of a button.

SEAL-RING OPTIONS

Use the Iceross Seal-In X Locking Liner with any of the three movable seal-ring options:



X-Classic
The classic
day-to-day seal



X-Grip
For good
rotational control



X-Volume
Accommodates more
volume fluctuation



Waterproof

Specifications

LOCK FEATURES:

- Option to change easily between locking suspension only, Seal-In only or combining these suspensions mechanisms at the same time.
- Simple release from liner with Quick Connect
- Waterproof
- Unity compatible
- Compatible with Direct Socket, 4-hole lamination sockets and laminated sockets



Icelock 850 Hybrid for Direct Socket



Icelock 850 Hybrid for 4-Hole Socket



Icelock 850 Hybrid for Laminated Socket

QUICK CONNECT FEATURES:

The Quick Connect is designed for ease and convenience, allowing users to don the socket while seated and release it from the liner with a simple push of a button.



Quick Connect w/ Icelock 850 Attachment Pin

ORDERING INFORMATION

Part#	Description
L-850001	Icelock 850 Hybrid for Direct Socket
L-850002	Icelock 850 Hybrid for 4-hole lamination
L-850003	Icelock 850 Hybrid for Laminated Sockets

The lock should be used with the Seal-In X Locking TF liner, Icelock 552 Expulsion Valve or Unity Expulsion valve, movable seal-rings and airtight socket.

Combine Seal-In Technology with Unity vacuum system for stabilization of volume fluctuations and comfortable socket fit throughout the day.

LINER FEATURES

- Locking and Seal-In suspension
- Silicone with active skin care ingredients
- Seamless functional textile
- Stabilising matrix
- Compatible with all movable seal-ring options
- Compatible with Unity



SEAL-IN X LOCKING



STABILIZING MATRIX



FABRIC COVER



ACTIVE SKIN CARE



SILKEN INNER SURFACE



CONICAL



EASY GLIDE

User Information	
Amputation Level:	Transfemoral
Impact Level:	Low to High
Liner Information	
Size Standard:	25,26.5,28,30,32,34,36,38,40,42,45,50,55
Size Conical:	23.5,25,26.5,28,30,32,34,36,38,40,42
Profile:	2mm Std and 3mm Con
Matrix:	Full Length
Size Seal Rings:	20,22,24,26,28,30,32,35,38,41,44,47,51,55, 60,65

PART NUMBER INFORMATION

Part#	Description	Contains
I-SXC0XX	Iceross Seal-In X-Classic	Seal-In X-Classic Ring
I-SXG0XX	Iceross Seal-In X-Grip	Seal-In X-Grip Ring
I-SXV0XX	Iceross Seal-In X-Volume	Seal-In X-Volume Ring
I-RXC001	Seal-In X-Classic Kit 26-51	10 Seal-In X-Classic Rings
I-RXC004	Seal-In X-Classic Kit 20-65	16 Seal-In X-Classic Rings
I-RXG001	Seal-In X-Grip Kit 26-51	10 Seal-In X-Grip Rings
I-RXG004	Seal-In X-Grip Kit 20-65	16 Seal-In X-Grip Rings
I-RXV001	Seal-In X-Volume Kit 26-51	10 Seal-In X-Volume Rings
I-RXV004	Seal-In X-Volume Kit 20-65	16 Seal-In X-Volume Rings