

ICEROSS® SYNERGY™ LOCKING LINER WITH WAVE




The Iceross Synergy liner with Wave incorporates a unique combination of two different silicone durometers to provide the perfect balance of comfort and stability.

A firm silicone outer layer provides outstanding stability for active use, while a softer shock-absorbing inner layer nurtures the skin, offering superb suspension and comfort.

Iceross Synergy with Wave makes the liner easier to flex at the knee, allowing it to adapt more to the user's movements. Less resistance means that less energy is required to flex the knee.

Össur recommends that Iceross Synergy Liner with Wave is used in conjunction with the Icelock® 600 Series.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Moderate to Extreme
LINER INFORMATION	
Size Standard:	16, 18, 20, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 40, 45.
Profile:	3mm
Matrix:	10cm or Custom 2-14cm*

SUSPENSION METHOD

FEATURES



LOCKING



MATRIX



FABRIC
COVER



ACTIVE SKIN
CARE



DOUBLE DURO



WAVE



MENTHOL



SILKEN INNER
SURFACE

PART NUMBER INFORMATION - ICEROSS SYNERGY LINER WITH WAVE

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-6303XX	3mm Wave	Locking	10cm
I-6206XX	3mm Wave	Locking	Custom 2-14cm*

*For a custom matrix, please specify the required matrix length on your order. Please note that deviation in the stabilizing matrix length can be +/- 1.5cm.


ICEROSS® SYNERGY™ LOCKING LINER



The Iceross Synergy liner incorporates a unique combination of two different silicone durometers to provide the perfect balance of comfort and stability.

A firm silicone outer layer provides outstanding stability for active use, while a softer shock-absorbing inner layer nurtures the skin, offering superb suspension and comfort.

Össur recommends that Iceross Synergy Liner is used in conjunction with the Icelock® 600 Series.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Moderate to Extreme
LINER INFORMATION	
Size Standard:	16,18,20,22,23.5,25,26.5,28,30,32,34,36,40,45.
Profile:	3mm
Matrix:	10cm or Custom 2-14cm*

SUSPENSION METHOD	FEATURES					
 LOCKING	 MATRIX	 FABRIC COVER	 ACTIVE SKIN CARE	 DOUBLE DURO	 MENTHOL	 SILKEN INNER SURFACE

PART NUMBER INFORMATION - ICEROSS SYNERGY LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-6003XX	3mm	Locking	10cm
I-6203XX	3mm	Locking	Custom 2-14cm*

*For a custom matrix, please specify the required matrix length on your order. Please note that deviation in the stabilizing matrix length can be +/- 1.5cm.