

# ICEROSS® SYNERGY™ CUSHION LINER WITH WAVE




Double the effect, with superb stability and comfort.

The Iceross Synergy Cushion liner with Wave incorporates a unique combination of two different silicone durometers to provide the perfect balance of comfort and stability.

A firm silicone outer layer provides outstanding stability for active use, while a softer shock-absorbing inner layer nurtures the skin, offering superb suspension and comfort.

The Wave feature makes the liner easier to flex at the knee, allowing it to adapt more to the user's movements. Less resistance means that less energy is required to flex the knee. It stretches easily over the patella and effectively controls folding of the liner behind the knee.

Össur recommends that Iceross Synergy Cushion Liner with Wave is used in conjunction with the Iceross Sleeve and the Icelock® 500 Series.



USER INFORMATION	
Amputation Level:	Transtibial
Impact Level:	Moderate to Extreme
LINER INFORMATION	
Size Standard:	16,18,20,22,23.5,25,26.5,28,30,32,34,36,40.
Profile:	3mm

## SUSPENSION METHOD



CUSHION

## FEATURES



FABRIC COVER



ACTIVE SKIN CARE



DOUBLE DURO



WAVE



MENTHOL



SILKEN INNER SURFACE

## PART NUMBER INFORMATION - ICEROSS SYNERGY CUSHION LINER WITH WAVE

PART#	PROFILE	SUSPENSION METHOD
I-CW63XX	3mm Wave	Cushion

# ICEROSS® SYNERGY™ CUSHION LINER




Double the effect, with superb stability and comfort.

The Iceross Synergy Cushion liner incorporates a unique combination of two different silicone durometers to provide the perfect balance of comfort and stability.

A firm silicone outer layer provides outstanding stability for active use, while a softer shock-absorbing inner layer nurtures the skin, offering superb suspension and comfort.

Össur recommends that Iceross Synergy Cushion Liner is used in conjunction with the Iceross Sleeve and the Icelock® 500 Series.



USER INFORMATION	
Amputation Level:	Transibial
Impact Level:	Moderate to Extreme
LINER INFORMATION	
Size Standard:	16,18,20,22,23.5,25,26.5,28,30,32,34,36,40.
Profile:	3mm

## SUSPENSION METHOD



CUSHION

## FEATURES



FABRIC COVER



ACTIVE SKIN CARE



DOUBLE DURO



MENTHOL



SILKEN INNER SURFACE

## PART NUMBER INFORMATION - ICEROSS SYNERGY CUSHION LINER

PART#	PROFILE	SUSPENSION METHOD
I-CL63XX	3mm	Cushion