


ICEROSS TRANSFEMORAL SEAL-IN® LINER



The Iceross Transfemoral Seal-In liner offers uncompromised security and freedom of movement for low to extremely active transfemoral amputees. Now with a new textile that has improved the adhesion by 25%, the fabric cover is extremely durable and allows radial stretch and comfortable elasticity.

Össur recommends that Iceross Seal-In Transfemoral Liner is used in conjunction with the Icelock® 500 Series Expulsion Valve.

USER INFORMATION	
Amputation Level:	Transfemoral
Impact Level:	All
LINER INFORMATION	
Size Standard:	25,26.5,28,30,32,34,36,38,40,45,50,55.
Size Conical:	25,26.5,28,30,32,34,36,38.
Profile:	3mm
Matrix:	Full Length



SUSPENSION METHOD



SEAL-IN

FEATURES



MATRIX



FABRIC COVER



ACTIVE SKIN CARE



SILKEN INNER SURFACE



CONICAL

PART NUMBER INFORMATION - ICEROSS TRANSFEMORAL SEAL-IN LINER

PART#	PROFILE	SUSPENSION METHOD	MATRIX
I-7532XX	Standard 3mm	Seal-In	Full Length
I-7632XX	Conical 3mm	Seal-In	Full Length