



ÖSSUR®

LIFE WITHOUT LIMITATIONS

MIAMI LUMBAR®

Because backs aren't flat.



MIAMI LUMBAR®.

THE BEST INNOVATE.

THE REST IMITATE.

Like an ergonomic chair.

The human anatomy is curved. Which is why both the Miami Lumbar belt itself and its various support panels are sculpted to hug the contours of the back, improving comfort, stabilization and compliance.



US

THEM




It's cool to be cool.

Breathable Aerospace fabric and ventilated support panels combine to keep Miami Lumbar cool and comfortable to wear for extended periods of time.

Pressure, only where you want it.

Miami Lumbar relieves pressure over the spinous processes while applying an even pressure to the paraspinal musculature to ensure comfort during healing.



It's a cinch to cinch.

With a four-to-one ratio drawstring and an easy-to-use pulley system, Miami Lumbar is simple to tighten, even for arthritic and/or geriatric patients.

A MODULAR SYSTEM OFFERING COMPLETE CONTROL

The trimmable Miami Lumbar® belt (slim and regular) combines with optional, heat-moldable anterior and posterior panels to create a customizable, modular system. step up, or down, the level of support to satisfy the immobilization needs of every patient.




The trimmable Miami Lumbar™ belt (slim shown above, regular shown below)



Posterior panels, with or without lateral supports, affix to the yellow attachment site (4 sizes)

Anterior panels (standard or pendulous) affix to the blue attachment site



Photos depict the highest level of immobilization: Miami Lumbar belt plus a pendulous anterior panel plus a posterior panel with lateral supports

TWO SIMPLE OPTIONS FOR ORDERING

OPTION 1: Order braces for inventory

UNIVERSAL MIAMI LUMBAR BELT		
		
SIZE	WAIST CIRC.	PART #
Univ Standard XS-2XL	26-50" (66-127cm)	LMB-006
Univ Large XS-6XL	26-66" (66-167cm)	LMB-010
Slim Univ Std XS-2XL	26-50" (66-127cm)	LMBSL-006
Slim Univ Large XS-6XL	26-66" (66-167cm)	LMBSL-010

ANTERIOR PANEL		
		
SIZE	WAIST CIRC.	PART #
Standard	<46" <116cm	LMB-ST
Slim	<46" <116cm	LMBSL-ST
Pendulous	>46" >116cm	LMB-PN

POSTERIOR PANEL WITHOUT LATERAL SUPPORTS		
		
SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	LMB-NS
Medium (M)	34-38" (86-96cm)	LMB-NM
Large (L)	38-50" (96-127cm)	LMB-NL
XLarge (XL)	50-66" (127-167cm)	LMB-NX

POSTERIOR PANEL WITH LATERAL SUPPORTS		
		
SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	LMB-LS
Medium (M)	34-38" (86-96cm)	LMB-LM
Large (L)	38-50" (96-127cm)	LMB-LL
XLarge (XL)	50-66" (127-167cm)	LMB-LX

OPTION 2: Order a brace for a specific patient

1. MIAMI LUMBAR BELT X X		
		
SIZE	WAIST CIRC.	XX CODE
XS	26-30" (66-76cm)	01
S	30-34" (76-86cm)	02
M	34-38" (86-96cm)	03
L	38-42" (96-106cm)	04
XL	42-46" (106-116cm)	05
2XL	46-50" (116-127cm)	06
3XL	50-54" (127-137cm)	07
4XL	54-58" (137-147cm)	08
5XL	58-62" (147-157cm)	09
6XL	62-66" (157-167cm)	10

2. POSTERIOR PANEL Y Y	
	
POSTERIOR PANEL DESC.	YY CODE
No Posterior Panel	NP
(1) Without Lateral Supports	NL
(2) With Lateral Supports	LS

3. ANTERIOR PANEL Z Z	
	
ANTERIOR PANEL DESC.	ZZ CODE
No Anterior Panel	NP
Standard Panel	ST
Pendulous Panel	PN

PART NUMBER = LMB-
FOR SLIM BELT ADD "SL" = LMBSL-0

XX

YY

ZZ

XX = Belt
YY = Posterior Panel
ZZ = Anterior Panel

NOTES

- Anterior and posterior panel sizes are automatically determined by belt size specified in Step #1
- Posterior panels are available in (S)mall, (M)edium, (L)arge and (XL)arge
- Anterior panels are available in standard (<46"/116cm) and pendulous (>46"/116cm)
- Adjustments to posterior panels may be necessary for waist circ. greater than 58"/147cm

FOLLOW ÖSSUR ON



USA
(800)233-6263
WWW.OSSUR.COM

CANADA
(800)663-5982
WWW.OSSUR.CA

