




MIAMI LUMBAR®



- 1 Both the Miami Lumbar belt and its various support panels are sculpted to hug the contours of the back, improving comfort, stabilization and compliance.
- 2 Breathable Aerospacer fabric and ventilated support panels combine to keep Miami Lumbar cool and comfortable.
- 3 Universal trimmable or sized belt options. Heat moldable anterior and posterior panels.
- 4 Miami Lumbar relieves pressure over the spinous processes while applying an even pressure to the paraspinal musculature to ensure comfort during healing.



SLIM MODEL ALSO AVAILABLE.



MIAMI LUMBAR® TLSO

- Building on the modularity of the Miami Lumbar LSO, the TLSO Kit helps you reduce your inventory, while enabling you to properly and confidently fit nearly every patient you see. The Miami Lumbar universal belt, a selection of panels, and a T-Connector Kit are all you need on-hand.
- Simply fit the Miami Lumbar LSO as you normally would. Then, set the desired height of the anterior and posterior thoracic extensions, adjust the twin-angle hinges and cinch down the trimmable shoulder straps to secure your patient in one of the most customizable off-the-shelf TLSO's available.
- As your patient progresses through different levels of healing, you can step-down the TLSO to a standard LSO brace, or remove the ATE and add the optional axillary shoulder straps to the PTE. Whatever level of immobilization you need to achieve, the Miami Lumbar TLSO accommodates you at every step.

PART#	ACCESSORIES
LMB-T01	TLSO system with ATE & PTE
LMB-T02	PTE and shoulder straps
LMB-A01	Air Bladder
LMB-A02	Pad replacement set

PART#	TLSO KITS
LMB-XXYYNPT	Standard Miami Lumbar TLSO Kit
LMBSL-0XXYYNPT	Slim Miami Lumbar TLSO Kit

XX = 01 (XS); 02 (S); 03 (M); 04 (L); 05 (XL);
06 (2XL); 07 (3XL); 08 (4XL); 09 (5XL); 10 (6XL)
YY = NL (Without Lateral Supports); LS (With Lateral Supports)

MIAMI LUMBAR® LSO & TLSO ORDERING INFORMATION

OPTION 1: To order braces for inventory

UNIVERSAL MIAMI LUMBAR BELT			ANTERIOR PANEL		
					
SIZE	WAIST CIRC.	PART #	SIZE	WAIST CIRC.	PART #
Univ Standard XS-2XL	26-50" (66-127cm)	LMB-006	Standard	<46" <116cm	LMB-ST
Univ Large XS-6XL	26-66" (66-167cm)	LMB-010	Slim	<46" <116cm	LMBSL-ST
Slim Univ Std XS-2XL	26-50" (66-127cm)	LMBSL-006	Pendulous	>46" >116cm	LMB-PN
Slim Univ Large XS-6XL	26-66" (66-167cm)	LMBSL-010	TLSO Kit	N/A	LMB-T01

POSTERIOR PANEL WITHOUT LATERAL SUPPORTS			POSTERIOR PANEL WITH LATERAL SUPPORTS		
					
SIZE	WAIST CIRC.	PART #	SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	LMB-NS	Small (S)	26-34" (66-86cm)	LMB-LS
Medium (M)	34-38" (86-96cm)	LMB-NM	Medium (M)	34-38" (86-96cm)	LMB-LM
Large (L)	38-50" (96-127cm)	LMB-NL	Large (L)	38-50" (96-127cm)	LMB-LL
XLarge (XL)	50-66" (127-167cm)	LMB-NX	XLarge (XL)	50-66" (127-167cm)	LMB-LX

OPTION 2: To order a brace for a specific patient

1. MIAMI LUMBAR BELT X X			2. POSTERIOR PANEL Y Y		3. ANTERIOR PANEL Z Z	
			 			
SIZE	WAIST CIRC.	XX CODE	POSTERIOR PANEL DESC.		ANTERIOR PANEL DESC.	
XS	26-30" (66-76cm)	01	No Posterior Panel		No Anterior Panel	
S	30-34" (76-86cm)	02	(1) Without Lateral Supports		Standard Panel	
M	34-38" (86-96cm)	03	(2) With Lateral Supports		Pendulous Panel	
L	38-42" (96-106cm)	04				
XL	42-46" (106-116cm)	05				
2XL	46-50" (116-127cm)	06				
3XL	50-54" (127-137cm)	07				
4XL	54-58" (137-147cm)	08				
5XL	58-62" (147-157cm)	09				
6XL	62-66" (157-167cm)	10				

PART NUMBER = LMB-    **FOR SLIM BELT ADD "SL-0" = LMBSL-0**

XX = Belt
YY = Posterior Panel
ZZ = Anterior Panel

NOTES

- Anterior and posterior panel sizes are automatically determined by belt size specified in Step #1
- Posterior panels are available in (S)mall, (M)edium, (L)arge and (XL)arge
- Anterior panels are available in standard (<46"/116cm) and pendulous (>46"/116cm)
- Adjustments to posterior panels may be necessary for waist circ. greater than 58"/147cm

1 Option not available when ordering TLSO Kit

2 When ordering TLSO Kit, use code LMB-XXYYNPT or LMBSL-0XXYYNPT.