NEXTSTEP

Flexible keel for the active K2 patient



- 2 degrees of Dorsiflexion during swing
- 15 degrees of Plantarflexion in stance
- Heel contact torque reduction



- Ensures knee stability
- **NEXTSTEP** promotes propulsion
- For the active K2 wearer

NEXTSTEP Keel Size / Category



- Category 4 Sizes 23 and 24 are rated to 80 kg or less than 175 lbs.
- Category 5 Sizes 25 and 26 are rated to 100 kg or less than 220 lbs.
- Category 6 Sizes 27 and 28 are rated to 125 kg or less than 275 lbs.

SIZE / CATEGORY	4	5	6
23 / 24			
25 / 26			
27 / 28			

NEXTSTEP with Ankle Glide:

- Comes with a foot shell, spectra sock
- Foot shell comes in two colors: beige & brown
- Ships standard with beige foot shell unless part number indicates brown foot shell.

When placing an order, please use the below formula.

1113-C (indicate category # see size/category chart)-size L or R For example, if you were ordering a size 23, left foot, category 4, the part number would be: 1113-C4-23L

NOTE: When ordering a **brown** foot shell, please add a "B" after the p/n: 1113-**B**-C4-23L

Foot Shell replacement

(BIEGE) = 1114-SH- Size(R/L)(BROWN) = 1114-B-SH-size(R/L)

PDAC APPROVED:

Assigned HCPCS Codes L5968 & L5972

FOOT INFORMATION		
Product #	1113	
Sizes:	23-28	
Combined weight of Foot Size 26	1.12 lbs (17.9 oz / 507.45g)	
Build Height to mid-dome of the ankle	3.8 inches. (96.5mm)	
Heel Height:	3/8" (10mm)	
Maximum Patient Weight	125kg / 275 lbs	
Warranty	2 year	

P/N	DESCRIPTION
4200	Foot Shell Removal Tool











Tel: (714) 524-0663 | Fax: (714) 364-8113

ST&G CORPORATION Rm. 911 43. Yangsan-ro

Yeongdeungpo-gu, Seoul, Korea

