

WillowWood One® TT System

Patient Criteria

The Patient criteria for using the WillowWood One System in a transtibial application are as follows:



Indications

Patients who meet the following criteria should be considered for this system:


- Patients who are transtibial amputees.
- Patients who comment about lack of suspension or feeling insecure in their current prosthesis.
- Patients who potentially would benefit from better limb health.
- Patients who exhibit high cognitive function and are “in tune” with prosthetic fit.
- Patients with reasonable hand dexterity.
- Community ambulators.
- Patients with redundant tissue.



Contraindications

Patients with any of the following conditions may not be considered for this system:

- Patients who are transfemoral amputees.
- Patients with a residual limb length less than 5 inches.
- Patients with bulbous residual limbs.
- Patients with open wounds.
- Patients who lack the ability to don and doff a gel interface or the ability to properly don and doff the system.
- Patients who have poor residual limb hygiene or who are unwilling or unable to comply with the cleaning and disinfecting instructions for the liner.
- Patients with a combination of limb length and limb circumference that is unsuitable according to the chart below.




WillowWood®

The Ohio Willow Wood Company
15441 Scioto Darby Road
Mt. Sterling, OH 43143
phone 740.869.3377 / 800.848.4930
fax 740.869.4374 www.willowwoodco.com



Ohio Willow Wood Company B.V.
Keizersgracht 62/64
1015 CS Amsterdam
The Netherlands



Patent www.willowwoodco.com/education-and-resources/patents

PN-2519-A 28 MARCH 2017