









A-Pylon

	Connection to modular components with pyramid adapter/receiver	Connection to modular components with tube clamp	Connection to socket adapter/attachment block
Accessories	 <p>4R82=P Tube Clamp Adapter Set Ø 34 mm SH* 12mm</p> <p>4R82 Tube Clamp Adapter Set Ø 34 mm SH* 33mm</p>	<p>Ø 30 mm Ø 34 mm</p>  <p>2R182=30 Transfemoral Fitting SH* 89 mm</p> <p>2R183 Spacer sleeve, 50 mm SH* 6 mm</p>	<p>4-hole Euro M6 4-hole Euro M6 without thread</p>  <p>4R431=1 Socket Adapter SH* 5 mm</p> <p>4R431=2 Socket Adapter SH* 5 mm</p>
	 <p>2R183 Spacer sleeve, 50 mm SH* 6 mm</p> <p>2R183=L Length adjustment, 120 mm SH* 79 mm</p>		 <p>4R415 Spacer Plate, 4-hole, height 3 mm SH* 3 mm</p>
Single components			

B-Pylon

	Connection to modular components with tube clamp	Connection to socket adapter/attachment block
Accessories	<p>Ø 30 mm Ø 34 mm</p>  <p>2R185=30 Transfemoral Fitting SH* 89 mm</p> <p>2R185=34 Transfemoral Fitting SH* 89 mm</p>	<p>4-hole Euro M6 4-hole Euro M6 without thread</p>  <p>4R432=1 Socket Adapter SH* 3 mm</p> <p>4R432=2 Socket Adapter SH* 3 mm</p>
Single components		 <p>4R415 Spacer Plate, 4-hole, height 3 mm SH* 3 mm</p>

* SH = system height