

Xtend® Foot Fri-Shock/Torsion H100 Technical Data Specification

The Xtend Foot Fri is a high-performance prosthetic foot with integrated shock absorption and torsional rotation, reducing impact and stress on the residual limb for smoother, more natural movement. It offers 17° of lateral flexibility for stability on uneven terrain, and its slim-profile design suits users with limited clearance. Made from highend adaptable glass and carbon fiber, it delivers excellent energy return and a low center of gravity - empowering confident, comfortable steps in both daily life and active environments.





Warranty

Xtend Foot Fri – 36 months. Footshell – 9 months.



Max weight patient

150 kg/330 lbs (size category P7).



One unit

Xtend Foot Fri has no loose parts.




Pre-assembled by delivery

Xtend Foot Fri, Footshell, Spectra sock already in place.

PR0027_US_03_2025-11-17 | We reserve the right to make technical changes.

Technical facts Xtend® Foot Fri H100





	Foot size (cm)	2	23	2	24	2	25	2	6	2	27	2	8	2	29
Category		А		А		Α		А	В	Α	В	А	В	А	
P3	Max 60 kg/130 lbs	72	142	70	142										
P4	Max 80 kg/175 lbs	72	142	70	143	76	151	82	153	84	154	86	155	87	156
P5	Max 100 kg/220 lbs	72	143	70	144	76	151	82	153	84	154	86	157	87	158
P6	Max 125 kg/275 lbs					76	153	82	153	84	154	86	159	87	159
D7	May 150 kg /220 lbs							0.2	152	0/	15.0	٥٤	160	07	1(1



Size guide

		Foot size (cm)									
Category		23	24	25	26	27	28	29			
P3	Max 60 kg/130 lbs	H100-323-STW	H100-324-STW								
P4	Max 80 kg/175 lbs	H100-423-STY	H100-424-STY	H100-425-STY	H100-426-STY	H100-427-STY	H100-428-STY	H100-429-STY			
P5	Max 100 kg/220 lbs	H100-523-STG	H100-524-STG	H100-525-STG	H100-526-STG	H100-527-STG	H100-528-STG	H100-529-STG			
P6	Max 125 kg/275 lbs			H100-625-STR	H100-626-STR	H100-627-STR	H100-628-STR	H100-629-STR			
P7	Max 150 kg/330 lbs				H100-726-STB*	H100-727-STB*	H100-728-STB*	H100-729-STB*			
* later launcl	٦										

How to order

For Caucasian, Brown, Tan or Jet Black footshell, add FC, FB, FT or FJ as a suffix after the article code. Example: H100-426-STY-FB.

For left or right footshell, add L or R as last part of the footshell color suffix. Example: H100-426-STY-FBL.

