



PTO Soft Knee Brace



Sizes: xs, s, m, l, xl, xxl

The PTO is designed to align the knee for patients with maltracking or patella tracking conditions. The PTO adjusts the amount of compression applied from the lateral pressure plate by forming a soft tissue wall that controls patellar subluxation.

Common Examples of Use

- Mild chondromalacia
- Patella subluxation
- Patella dislocation
- Patella tendonitis
- Chronic patellar maltracking
- Patellofemoral pain syndrome
- Post lateral release
- Post patellar realignment

Product Features

- Available in Airmesh or neoprene
- Open back or closed back options available
- Available with aluminum high performance hinge bars
- Dynamic Function for normal gait

Category: Patellofemoral.

Instructions for Use

Instructions for Use

Part Numbers

Part Numbers

Part #	Description
1419X	PTO, Airmesh, Left
1418X	PTO, Airmesh, Right
1423X	PTO, Airmesh, Open Back, Left
1422X	PTO, Airmesh, Open Back
1431X	PTO, Airmesh, High Performance, Open Back, Left
1430X	PTO, Airmesh, High Performance, Open Back, Right
1417X	PTO, Neoprene, Left
1416X	PTO, Neoprene, Right
1421X	PTO, Neoprene, Open Back, Left
1420X	PTO, Neoprene, Open Back, Right
7181X	PTO High Performance Sports Cover
7179X	PTO HP Replacement Bar Kit, Left
7180X	PTO HP Replacement Bar Kit, Right

Product Sizing Chart

Image not found or type unknown

