



RANGE OF MOTION (ROM) WALKER

Indications: Stable foot and ankle fractures, severe sprains or strains, soft tissue injuries, post-op and post-cast removal

- Range of Motion (ROM) joint can be fixed in any position between 0° - 45° of plantar flexion
- The ROM joint adjusts in 7.5° increments in both dorsi and plantar flexion and can either be locked out at any desired interval or left with desired range
- Cushioned insole and heel for increased comfort
- Malleable aluminum uprights for a customized patient fit
- Suggested HCPCS Code: L4386 or L4387

Size	High Top	Low Top	Men	Women
Small	T14022	T14122	4-6	5-7
Medium	T14024	T14124	7-10	8-11
Large	T14026	T14126	11-13	12-14



PEDIATRIC FRACTURE WALKER

Indications: For stable ankle sprains and fractures, forefoot sprains and stable distal tibial fractures, and soft tissue injuries

- Steel enforced plastic molded uprights are slotted for easier and more secure strapping
- Heel height, width and rocker bottom sole designed for small children
- Shock-absorbing insole reduces impact of heel strike when walking
- Suggested HCPCS Code: L4387

Size	Age	Foot Length	Part #
Small/Medium	1 - 2 years	Up to 6"	T13022
Large	2.5 - 6 years	6" - 7.5"	T13024
X-Large	6.5 - 9.5 years	7.5" - 8.5"	T13026