

PRO SELECT WALKER



Indications: soft tissue injury; grades 2 & 3 strains and sprains, stable fractures

- Ultra low heel height in conjunction with the rocker bottom sole encourages a more natural gait pattern and reduces pelvic tilt which mitigates secondary hip and back pain
- Patent pending technology allows air to circulate to prevent the accumulation of heat during patient wear
- Small perforations work to prevent humidity and help heat escape instead of accumulate under the foot
- Flared plastic uprights help reduce ankle rotation
- Extra wide foot bed to accommodate for patient dressings
- Suggested HCPCS Code: L4387



Size	High Top	Low Top	Men	Women
Small	T17022	T17122	4.5-7	6-8
Medium	T17024	T17124	7.5-10	8.5-11.5
Large	T17026	T17126	10.5-12.5	12-13.5
X-Large	T17028	T17128	13+	14+

*Liner and strap extenders available by request

PRO SELECT PNEUMATIC WALKER

Indications: soft tissue injury; grades 2 & 3 strains and sprains, stable fractures.

- Ultra low heel height in conjunction with the rocker bottom sole encourages a more natural gait pattern and reduces pelvic tilt which mitigates secondary hip and back pain
- Patent pending technology allows air to circulate to prevent the accumulation of heat during patient wear
- Small perforations work to prevent humidity and help heat escape instead of accumulate under the foot
- Flared plastic uprights help reduce ankle rotation
- Reinforced softgood liner to prevent disintegration during use
- Easy to use, extra durable turn knob system to inflate and deflate the pneumatic liner
- Suggested HCPCS Code: L4361



Size	High Top	Low Top	Men	Women
X-Small	T17200	T17300	<4	<5.5
Small	T17222	T17322	4.5-7	6-8
Medium	T17224	T17324	7.5-10	8.5-11.5
Large	T17226	T17326	10.5-12.5	12-13.5
X-Large	T17228	T17328	13+	14+

*Liner and strap extenders available by request

