

Conversion table	
Kg	Lbs
45	99
80	176
100	220
125	275
150	330



• Pylon without Cup

To be adapted on a cup. Can be assembled to the distal end of tube 1G06.

Item N°.	Max. Patient's Weight	Material	Weight	TH	Diameter
1G2648	100 kg	Aluminium	225 g	480 mm	25 mm



• Rubber Cup for Pylon

Ø 30 mm base.

Item N°.	Max. Patient's Weight	Material	Weight	TH	Connection
EE022	100 kg	Rubber	140 g	65 mm	Ø 25 mm 1G2648



• Adapter Ring

Made for fitting components on tubes of different diameters.

Item N°.	Max. Patient's Weight	Material	Weight	TH	Diameter
1X68	100 kg	Aluminium	18 g	26 mm	from 25 to 22 mm
1X69	100 kg	Aluminium	29 g	34 mm	from 30 to 22 mm
1X70	100 kg	Aluminium	20 g	34 mm	from 30 to 25 mm
1X71	100 kg	Aluminium	22 g	39 mm	from 34 to 30 mm



• Tube Ø 30 mm with Pyramid Receiver

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1D52	100 kg	Aluminium	111 g	130 mm	63 mm



• Tube Ø 34 mm with Pyramid Receiver

Item N°.	Max. Patient's Weight	Material	Weight	TH	BH
1D52-P6	125 kg	Aluminium	121 g	130 mm	67 mm