



Many patients lack the room for a distal lock, feel uncomfortable in a pull-in suction socket or have difficulty engaging a distal pin, particularly older patients who require assistance donning their prosthesis. In these instances, the Summit™ Lock is ideal for patients of any activity level or body weight. It turns a cushion liner into a locking liner with no clearance required. Controls rotation, easy to use even on existing sockets and can be glued to most fabric liners.

#### Summit™ Lock Kit

Order No.	Description
<b>CD101</b>	Standard
<b>CD102</b>	Low Profile

#### Summit™ Lock Parts and Accessories

Order No.	Part Description
<b>CD102XT</b>	Extra Tab Kit ( <i>for additional liner</i> )
<b>CD4120</b>	Coyote Quick Adhesive, 20cc
<b>CD4747</b>	Coyote Quick Adhesive, 50ml
<b>CDDMA51</b>	Dispensing Gun for Adhesive
<b>CD401</b>	Mix Tips for Adhesive, 6 Pack
<b>CD402</b>	Mix Tips for Adhesive, 20 Pack