

A full-page photograph of a middle-aged man running on a paved path. He is wearing a bright yellow-green jacket over a blue t-shirt, blue shorts, a black cap, and black running shoes. He has a white and blue Bauerfeind knee brace on his left knee. The background shows a large tree and a clear blue sky. Two horizontal bars, one blue and one yellow, are positioned above and below the main text respectively.

# KNEE OA SOLUTIONS – SUPPORTING AN ACTIVE AND PAIN-FREE LIFESTYLE

# KNEE OSTEOARTHRITIS

## A MAJOR CHALLENGE IN ORTHOPEDICS.

Degenerative arthritis of the knee (osteoarthritis) is a condition in which the natural cushion in between the joints “wears” away. According to the Arthritis Foundation more than 27 million Americans have osteoarthritis, which can occur in almost any joint but is most common in the knee.

There are a range of treatment options available to help manage symptoms associated with knee OA. Common modalities include physical therapy, weight loss, and activity modification. Injections (visco-supplementation) are commonly used for more severe knee OA.

Knee braces and other supportive devices such as cushioning gel insoles can often help relieve pain and provide joint stability.

Finally, if osteoarthritis symptoms are severe and other treatments don’t succeed, knee replacement surgery may be required.

### **PAIN RELIEF**

HIGH LEVEL OF WEARING COMFORT FOR -

### **BEST PATIENT COMPLIANCE**

EASY TO FIT WITH -

### **COST EFFECTIVE REIMBURSEMENT**

# SecuTec® OA

## ORTHOSIS FOR UNLOADING AND STABILIZATION OF THE KNEE JOINT.

### INDICATIONS

- Moderate to severe medial knee osteoarthritis (OA)
- Moderate to severe lateral knee osteoarthritis (OA)
- Post-operative stabilization and unloading (relief)

### TOOL-FREE:

- Extension 0°, 10°, 20°, 30°, 45°
- Flexion 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Fixation 0°, 10°, 20°, 30°, 45°  
(optional – not included with delivery)

**PDAC verified**  
**L1845/L1852**

**NEW**

Anatomically  
contoured  
condylar  
pads

Tool-free  
hinge  
adjustments

Medial and lateral  
hinges for optimum  
varus and valgus  
adjustment



Distal frame adjustments allow  
for optimum force transmission  
without compromising knee width





# BEST WEARING COMFORT

## EASY TO PUT ON. COMFORTABLE. LIGHTWEIGHT.



### OPTIMUM PRESSURE DISTRIBUTION

#### 01

- Extra-wide condylar pads for added comfort
- Cushioning wedges to ensure a flat fit

#### 02

- Extra side pad relieves the head of the fibula



### EXCELLENT WEARING COMFORT

- Anti-shock technology provides a more fluid walking movement
- Microfiber cushioning
- High-grade aluminum alloy makes the frame extremely lightweight
- Quick snap closure system for convenience and consistent tension



### PERFECT FIT

- Dynamic cuff adjusts to the calf muscle for optimum pressure distribution

# BAUERFEIND KNEE OA SOLUTIONS FOR ALL WALKS OF LIFE.

Bauerfeind is known globally for world class quality, innovation and manufacturing of orthopedic braces that provide proven therapeutic benefits and are extremely comfortable to wear. With the new SecuTec OA, Bauerfeind now offers a complete range of knee OA solutions covering the different stages of OA progression.

			 <div>NEW</div>	
Classification based on Kellgren-Lawrence scale	GRADE 1	GRADE 2	GRADE 3	GRADE 4
<b>GenuTrain</b>	<b>GenuTrain A3</b>	<b>GenuTrain S/S Pro L1832/L1833</b>	<b>SecuTec OA L1845/L1852</b>	
Mild osteoarthritis	Mild to moderate osteoarthritis of the knee without leg misalignment	Mild to moderate osteoarthritis of the knee with joint insta- bility	Moderate to severe osteoarthritis of the knee with varus or valgus misalignment	

## Get this brace and get off the bench.

This product is genuine. Don't let the price point turn you away. It's worth every penny. I have a serious case of osteoarthritis in my knee and this brace gave me tremendous support, relief and mobility. If you're looking for an easy-to-use brace that will give you, your active lifestyle back, this is the one ... leave all others behind.

Online review, patient using GenuTrain S

**BAUERFEIND USA, INC.**

75 14th St NE

Suite 2350

Atlanta, GA 30309

P 1 800-423-3405

P 404-201-7800

F 404-201-7839

[info@bauerfeindusa.com](mailto:info@bauerfeindusa.com)



Bauerfeind supports are used and endorsed  
by the National Basketball Athletic Trainers  
Association.