


SILITOUCH Aloe Locking TT			
Item N°.	Thickness	Distal thickness	Size
1S711-3-160	3 mm	17.5 mm	16
1S711-3-180			18
1S711-3-200			20
1S711-3-220			22
1S711-3-235			23.5
1S711-3-250			25
1S711-3-265			26.5
1S711-3-280			28
1S711-3-300			30
1S711-3-320			32
1S711-3-340			34
1S711-3-360			36
1S711-3-400		18.5 mm	40
1S711-3-450		18.5 mm	45
1S711-6-160	6 mm	19 mm	16
1S711-6-180			18
1S711-6-200			20
1S711-6-220			22
1S711-6-235			23.5
1S711-6-250			25
1S711-6-265			26.5
1S711-6-280			28
1S711-6-300			30
1S711-6-320			32
1S711-6-340			34
1S711-6-360			36
1S711-6-400		20 mm	40
1S711-6-450		20 mm	45

SILITOUCH Aloe Cushion TT			
Item N°.	Thickness	Distal thickness	Size
1S721-3-160	3 mm	14.5 mm	16
1S721-3-180			18
1S721-3-200			20
1S721-3-220			22
1S721-3-235			23.5
1S721-3-250			25
1S721-3-265			26.5
1S721-3-280			28
1S721-3-300			30
1S721-3-320			32
1S721-3-340			34
1S721-3-360			36
1S721-3-400		16 mm	40
1S721-3-450		16 mm	45
1S721-6-160	6 mm	16 mm	16
1S721-6-180			18
1S721-6-200			20
1S721-6-220			22
1S721-6-235			23.5
1S721-6-250			25
1S721-6-265			26.5
1S721-6-280			28
1S721-6-300			30
1S721-6-320			32
1S721-6-340			34
1S721-6-360			36
1S721-6-400		17.5 mm	40
1S721-6-450		17.5 mm	45

SILITOUCH Classic Locking TT			
Item N°.	Thickness	Distal thickness	Size
1S712-3-160	3 mm	17.5 mm	16
1S712-3-180			18
1S712-3-200			20
1S712-3-220			22
1S712-3-235			23.5
1S712-3-250			25
1S712-3-265			26.5
1S712-3-280			28
1S712-3-300			30
1S712-3-320			32
1S712-3-340			34
1S712-3-360			36
1S712-3-400		18.5 mm	40
1S712-3-450		18.5 mm	45
1S712-6-160	6 mm	19 mm	16
1S712-6-180			18
1S712-6-200			20
1S712-6-220			22
1S712-6-235			23.5
1S712-6-250			25
1S712-6-265			26.5
1S712-6-280			28
1S712-6-300			30
1S712-6-320			32
1S712-6-340			34
1S712-6-360			36
1S712-6-400		20 mm	40
1S712-6-450		20 mm	45

SILITOUCH Classic Cushion TT			
Item N°.	Thickness	Distal thickness	Size
1S722-3-160	3 mm	14.5 mm	16
1S722-3-180			18
1S722-3-200			20
1S722-3-220			22
1S722-3-235			23.5
1S722-3-250			25
1S722-3-265			26.5
1S722-3-280			28
1S722-3-300			30
1S722-3-320			32
1S722-3-340			34
1S722-3-360			36
1S722-3-400		16 mm	40
1S722-3-450		16 mm	45
1S722-6-160	6 mm	16 mm	16
1S722-6-180			18
1S722-6-200			20
1S722-6-220			22
1S722-6-235			23.5
1S722-6-250			25
1S722-6-265			26.5
1S722-6-280			28
1S722-6-300			30
1S722-6-320			32
1S722-6-340			34
1S722-6-360			36
1S722-6-400		17.5 mm	40
1S722-6-450		17.5 mm	45



HOW TO DETERMINE THE SIZE  
Metric measurement of the residual limb circumference, at 4 cm from the distal end.

THE PROTEOR PROMISE

PROTEOR is a leading prosthetic and orthotic device company committed to collaborating with clinicians to spark real-world innovations. We are inspired by the people who use our products. Seeing their everyday achievements motivate us to make more things possible. Our extensive portfolio includes technological breakthroughs that allow the utmost mobility, maneuverability, adaptability, and durability. PROTEOR's process for developing advancements and community support is based on a collaborative approach with clinicians and the people they care for. When we all work together, it's amazing how much patients can achieve.

CONNECT WITH US  
f i y t i n


US Headquarters | 1236 West Southern Ave, Suite 101,  
Tempe, AZ 85282 855.450.7300 us.proteor.com



SILITOUCH

LINERS FOR LIVING

DESIGNED TO DELIVER ALL-DAY COMFORT





# SILITOUCH

## RELIABLE COMFORT.

Counting on comfort, delivering durability



### A NEW DIRECTION IN COMFORT

PROTEOR engineers took patient complaints about the damage of textile tears to heart. By strengthening the fabric and adding a distal matrix, their approach means finishing the day without worry of liner failures.



### COMFORT THAT SOOTHES

Caring for your patient's skin is critical to preventing activity-ending sores. With an Aloe option, you can add extra comfort to complement the durable outer cover.



### PROTECTION THAT LASTS

Gravity dictates that the distal end of the residual limb endures force with every step. Based on the liner size, the silicone layer increases up to 20mm to counter everyday action.



### STYLE THAT SHINES

The elegant look and smooth touch of the liner fabric make it a pleasure to wear (and show off.) Easy to care for, simply hand-wash and dry and it's ready for another day.



### KEEPING ON

The fabric is placed to reduce tears that run, increasing usable life and reducing day-ending disappointment.

### MATRIX MATTERS

The locking liner provides reliable suspension thanks to its enhanced outer fabric with a 10cm distal matrix.




### 3MM OR 6MM


Your choice of liner thickness.

### CUSHION THAT COUNTS

From 14mm to 20mm, the distal silicone thickness varies based on the size and style, designed to match individual needs.



Option with distal attachment



Option without distal attachment

“ I JUST PUT IT ON AND GO. STARTING AT 6:30 AM TO 10 PM, I'VE NOTICED THAT THERE IS NO SLIPPING. I'LL WEAR THIS ONE INTO THE GROUND.”

LANCE K.