

PRACTITIONER SERIES

Defining Roles in Your stWear® Account Dashboard



Spinal.Tech/stWear

This e-book is designed as a companion piece to our Online Video Training.




Visit [YouTube.com/@SpinalTechnology](https://www.youtube.com/@SpinalTechnology) to watch the full video.

Primary Admin

Each account has one designated Primary Admin who oversees all Users*, physicians, and patients. This User:

- Has the ability to add or remove additional Admin Users within the account.
- Provides a single, unchangeable email address for the account upon setup.
- Can add or remove Users, including Admins, Practitioners, and Physicians.
- Can also be set up as a Practitioner.
- Has the ability to add or delete patients and assign devices.
- Can log in to the stWear® app to take readings.
- Cannot be set up as a Practitioner or Physician—this is strictly an administrative role.

*A User is defined as another Admin, Practitioner, or Physician.



← Users

Add new user

Cancel Save

New user details

First name Travis ✓


Last name Hood ✓

Email travis@spinaltech.com ✓

Admin ☒

Practitioner ☐

Physician ☐



Email travis@spinaltech.com ✓

Admin ☒

Practitioner ☒


Physician ☐

Practitioner

With the Practitioner role*, a User can:

- Hold both Practitioner and Admin roles (recommended for full account access).
- Add or delete only their assigned patients.
- See all patients within the account.
- Add or remove a patient's referring Physician.

*A User cannot be designated as a Physician if already set up as a Practitioner.



Patients

Users - ST Clinic Test

+ Add user


Search in all columns...			
Name	Email	Permission level	Status
Doctor Smith	doctorsmith@spinaltech.com	Physician	Active
Jake Mele	jake_me@spinaltech.com	Practitioner Admin	Active
Jenny Netz	jenny@stbauer.com	Practitioner	Active
Clinician Test	physician2@spinaltech.com	Admin Practitioner	Invited
Test Physician	physician@spinaltech.com	Physician	Active
Practitioner Test	practitioner@spinaltech.com	Primary Admin	Active

Physician

With the Physician role*, a User can:

- Only view data for patients assigned to them.
- View patient data within the dashboard only—the mobile app is not available for Physicians.

*This User does not have the ability to add or remove patients, assign devices, or modify account settings.



Name	Email	Permission level	Status
	practitioner@spinaltech.com	Primary Admin	Active
Travis Hood	travis@spinaltech.com	Practitioner Admin	Invited

If you need assistance, visit [Spinal.Tech/stWear](https://spinal.tech/stWear) to submit a support ticket or email stWear@SpinalTech.com.



Spinal
Technology

Spinal.Tech/**stWear**

Phone: (508) 957-8282

