

Carbon-prepreg Ankle-Foot Orthosis with enhanced calf stabilization, offering dynamic support and includes two sets of cushioning pads for added comfort.



SIZE GUIDE

Size	Side	Shoe size	Height	Item-No.
S	left	US Women's 5.5 – 8	34 cm	197R6/L1
M	left	US Women's 8 – 10.5 US Men's 6.5 – 9	35 cm	197R6/L2
L	left	US Men's 9 – 11.5	35 cm	197R6/L3
XL	left	US Men's 11.5 – 14	36 cm	197R6/L4
-	-		-	-
S	right	US Women's 5.5 – 8	34 cm	197R6/R1
M	right	US Women's 8 – 10.5 US Men's 6.5 – 9	35 cm	197R6/R2
L	right	US Men's 9 – 11.5	35 cm	197R6/R3
XL	right	US Men's 11.5 – 14	36 cm	197R6/R4

The PeroSupport.pro may be customized to foot shape and size by cutting the sole as per requirements. The supplied micro-fleece insole can be cut into shape as well. Carbon-prepreg is not thermo-formable.

The positive effect is only guaranteed if the orthosis is worn together with a shoe.

INTENDED USE

- ▶ Dynamic readjustment of the foot
- ▶ Restriction of plantar flexion and support of dorsal flexion
- ▶ Prevents lateral ankle twisting (supination position)
- ▶ Prevents side-slips (lateralization)
- ▶ Great energy storage and return by carbon-prepreg material

INDICATIONS

- ▶ Peroneal paresis
- ▶ Weakness of dorsiflexion of the foot
- ▶ Malalignment of the ankle joint as a result of weakness of dorsiflexion of the foot

CONTRAINDICATIONS

- ▶ Material interactions/allergies
- ▶ Do not wear directly on injured, irritated, or damaged skin
- ▶ Weight limit: 120 kg
- ▶ Deformation of the foot or ankle
- ▶ Swelling and extension of the calf skin
- ▶ Sensory insensitivity (e.g., diabetic polyneuropathy)
- ▶ Diabetic foot syndrome/ulcer
- ▶ Severe circulatory disorders/congestive edema

