



STUMP SHRINKERS

Comfortable medi knit technology

COMFORTABLE:

- Thermoregulating yarn with moisture wicking helps to keep limb cool and dry
- Antibacterial properties minimize odor
- Smooth distal seam
- Smooth transition from stocking to TF hip attachment

RELIABLE:


- Integrated silicone beaded top band prevents TT roll down
- Customized hip attachment positioning allows precise TF closure placement

EFFECTIVE:

- Two levels of compression: 20-30 mmHg and 30-40 mmHg
- Gradient compression from distal to proximal helps channel excess fluids out of limb

STUMP SHRINKERS

TRANSFEMORAL:




¹ at the thickest area of the thigh
² 5 cm above the distal end of the residual limb

Length	Size	Circumference cF²	Circumference cG¹	Part No.
Short 18 – 22 cm	1	30 – 32 cm	42 – 45 cm	AKXS0W1
	2	32 – 35 cm	45 – 49 cm	AKXS0W2
	3	35 – 38 cm	49 – 53 cm	AKXS0W3
	4	38 – 41 cm	53 – 57 cm	AKXS0W4
	5	41 – 44 cm	57 – 61 cm	AKXS0W5
	6	44 – 47 cm	61 – 65 cm	AKXS0W6
Medium 23 – 27 cm	1	30 – 32 cm	42 – 45 cm	AKXM0W1
	2	32 – 35 cm	45 – 49 cm	AKXM0W2
	3	35 – 38 cm	49 – 53 cm	AKXM0W3
	4	38 – 41 cm	53 – 57 cm	AKXM0W4
	5	41 – 44 cm	57 – 61 cm	AKXM0W5
	6	44 – 47 cm	61 – 65 cm	AKXM0W6
Long 28 – 32 cm	1	30 – 32 cm	42 – 45 cm	AKXL0W1
	2	32 – 35 cm	45 – 49 cm	AKXL0W2
	3	35 – 38 cm	49 – 53 cm	AKXL0W3
	4	38 – 41 cm	53 – 57 cm	AKXL0W4
	5	41 – 44 cm	57 – 61 cm	AKXL0W5
	6	44 – 47 cm	61 – 65 cm	AKXL0W6
X-Long 33 – 37 cm	1	30 – 32 cm	42 – 45 cm	AKXX0W1
	2	32 – 35 cm	45 – 49 cm	AKXX0W2
	3	35 – 38 cm	49 – 53 cm	AKXX0W3
	4	38 – 41 cm	53 – 57 cm	AKXX0W4
	5	41 – 44 cm	57 – 61 cm	AKXX0W5
	6	44 – 47 cm	61 – 65 cm	AKXX0W6

Select Compression Class to complete Part Number: X = 1 = 20 – 30 mmHg, X = 2 = 30 – 40 mmHg

TRANSTIBIAL:




³ at the middle of the thigh
⁴ at mid-patella
⁵ 5 cm above the distal end

Length	Size	Circumference cC⁵	Circumference cE⁴	Circumference cF³	Part No.
Short 28 – 34 cm	0	26 – 28 cm	28 – 30 cm	38 – 40 cm	BKXS0N0
	1	28 – 30 cm	30 – 32 cm	40 – 42 cm	BKXS0N1
	2	30 – 33 cm	32 – 35 cm	42 – 45 cm	BKXS0N2
	3	33 – 36 cm	35 – 38 cm	45 – 48 cm	BKXS0N3
	4	36 – 39 cm	38 – 41 cm	48 – 51 cm	BKXS0N4
	5	39 – 42 cm	41 – 44 cm	51 – 54 cm	BKXS0N5
Medium 35 – 41 cm	0	26 – 28 cm	28 – 30 cm	38 – 40 cm	BKXM0N0
	1	28 – 30 cm	30 – 32 cm	40 – 42 cm	BKXM0N1
	2	30 – 33 cm	32 – 35 cm	42 – 45 cm	BKXM0N2
	3	33 – 36 cm	35 – 38 cm	45 – 48 cm	BKXM0N3
	4	36 – 39 cm	38 – 41 cm	48 – 51 cm	BKXM0N4
	5	39 – 42 cm	41 – 44 cm	51 – 54 cm	BKXM0N5
Long 42 – 48 cm	0	26 – 28 cm	28 – 30 cm	38 – 40 cm	BKXL0N0
	1	28 – 30 cm	30 – 32 cm	40 – 42 cm	BKXL0N1
	2	30 – 33 cm	32 – 35 cm	42 – 45 cm	BKXL0N2
	3	33 – 36 cm	35 – 38 cm	45 – 48 cm	BKXL0N3
	4	36 – 39 cm	38 – 41 cm	48 – 51 cm	BKXL0N4
	5	39 – 42 cm	41 – 44 cm	51 – 54 cm	BKXL0N5
X-Long 47 – 53 cm	0	26 – 28 cm	28 – 30 cm	38 – 40 cm	BKXX0N0
	1	28 – 30 cm	30 – 32 cm	40 – 42 cm	BKXX0N1
	2	30 – 33 cm	32 – 35 cm	42 – 45 cm	BKXX0N2
	3	33 – 36 cm	35 – 38 cm	45 – 48 cm	BKXX0N3
	4	36 – 39 cm	38 – 41 cm	48 – 51 cm	BKXX0N4
	5	39 – 42 cm	41 – 44 cm	51 – 54 cm	BKXX0N5
XX-Long 52 – 58 cm	0	26 – 28 cm	28 – 30 cm	38 – 40 cm	BKXX0N0
	1	28 – 30 cm	30 – 32 cm	40 – 42 cm	BKXX0N1
	2	30 – 33 cm	32 – 35 cm	42 – 45 cm	BKXX0N2
	3	33 – 36 cm	35 – 38 cm	45 – 48 cm	BKXX0N3
	4	36 – 39 cm	38 – 41 cm	48 – 51 cm	BKXX0N4
	5	39 – 42 cm	41 – 44 cm	51 – 54 cm	BKXX0N5

Select Compression Class to complete Part Number: X = 1 = 20 – 30 mmHg, X = 2 = 30 – 40 mmHg

KEY FEATURES:



CUSTOMIZED HIP POSITIONING

Precise closure placement for transfemoral care



COMFORTABLE TO WEAR

Smooth transition from stocking to hip attachment



HIGH COMFORT

Smooth distal seam



RELIABLE GRIP

Silicone beaded top band for transtibial care

FOLLOW ÖSSUR ON



TEL (800) 233-6263
CANADA (800) 663-5982
WWW.OSSUR.COM

ÖSSUR
LIFE WITHOUT LIMITATIONS