

## USA Walker – Standard and Ankle

The USA Walker provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Fixed, aluminum uprights allow for the frame to be adjusted to conform better to the contour of the leg.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.

Part #	Description		
1200X	USA Walker – Standard		
1201X	USA Walker – Ankle		
1200XR	USA Walker Replacement Liner – Standard		
1201XR	USA Walker Replacement Liner – Ankle		
WLE	Walker Liner Enlarger – One Size		

  

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
3	SM	5 – 6	5.5 – 7
5	MD	6.5 – 10	7.5 – 11
7	LG	10.5 – 12	11.5+
9	XL	12.5+	–



## USA Air Walker – Standard and Ankle

The USA Air Walker provides maximum stability and support following injury or post-operative procedures.

### EXAMPLES OF USE

- Acute ankle sprains
- Soft tissue injuries
- Stable lower leg, foot and ankle fractures
- Post-operative use

### PRODUCT FEATURES

- Fixed, aluminum uprights allow for the frame to be adjusted to conform better to the contour of the leg.
- Medial/lateral air bladders offer customized compression to increase comfort and reduce edema.
- Shock-absorbing insole reduces impact of heel strike when walking.
- Rocker sole helps promote a natural gait to allow continuation of daily activities.



Part #	Description		
1210X	USA Air Walker – Standard		
1211X	USA Air Walker – Ankle		
L1210X	USA Air Walker Replacement Liner – Standard		
NSL1211X	USA Air Walker Replacement Liner – Ankle		
WLE	Walker Liner Enlarger – One Size		

X =	Size	Men's Shoe Size	Women's Shoe Size
1	XS	–	3 – 5
3	SM	5 – 6	5.5 – 7
5	MD	6.5 – 10	7.5 – 11
7	LG	10.5 – 12	11.5+
9	XL	12.5+	–