

Wrist-Hand-Thumb Orthosis Anatomically shaped of low density polyethylene with closed-cell foam lining. Tongue is rigid yet comfortable for most traumatic, post-surgical, or longer term rehabilitation applications. Adjustable Velcro® straps allow for swelling and facilitate wound management. Easily trimmed with cast scissors. May be reshaped with heat gun.

- 2 socks (with thumb)
- · Large thumb component

Padded Tongue

· Perforated for ventilation

Model #

Size	Left	Right	Palm Width	Forearm Circ.	Brace Length
Small	1000	1005	21/2"-3"	7"-9"	71/2"
Medium	1001	1006	3"-1/2"	9"-101/2"	7¾"
Large	1002	1007	3½"-4"	10"-12"	81/2"

Available as: Lined WHTO w/o socks (.02 series) e.g., 1001.02

Standard model w/o socks or lining (.01 series) e.g., 1001.01

Replacement socks (2) (Model # 4135)