

28U33

## ***WalkOn Lateral.***

The *WalkOn Lateral* lifts the foot and also stabilizes the ankle joint.

### **Indication**

Drop foot, with or without mild spasticity,  
e.g. caused by:


- Peroneal paralysis
- Stroke
- Traumatic brain injury
- Multiple sclerosis
- Neuromuscular atrophy

The indication must be determined by the doctor.

### **Benefits/advantages**

- High energy return of the orthosis
- Effectively supports a physiological gait, even at different speeds and with different loads
- Stabilizes the ankle joint
- Increases the stride length by improving muscle dynamics
- The orthosis allows weight to be placed on the forefoot for climbing stairs and squatting
- Limits plantar flexion and supination when worn with a sturdy shoe
- Slim, lightweight design
- Optimized compliance thanks to a removable and washable pad/strap combination that can be aligned in the preferred direction for tightening the strap
- Optimized durability: adult sizes tested for two million steps and junior sizes tested for one million steps
- Material: carbon fiber prepreg
- Available in a variety of adult and junior sizes

**ottobock.**





### 28U33 WalkOn Lateral

Article number	Side	Shoe size (US)	Calf height
28U33=L36-39	Left	Men: 5-6	12.8" (32.5 cm)
28U33=R36-39	Right	Women: 6-8.5	
28U33=L39-42	Left	Men: 6-9	13.8" (35 cm)
28U33=R39-42	Right	Women: 8.5-10	
28U33=L42-45	Left	Men: 9-12	14.6" (37 cm)
28U33=R42-45	Right	Women: 10-13	
28U33=L45-48	Left	Men: 12-15	15.4" (39 cm)
28U33=R45-48	Right		

### 28U33 WalkOn Lateral Junior

Article number	Side	Shoe size (US)	Calf height
28U33=L24-27	Left	Junior: 7.5/8-10	8" (20.5 cm)
28U33=R24-27	Right		
28U33=L27-30	Left	Junior: 10-12.5	9.3" (24 cm)
28U33=R27-30	Right		
28U33=L30-33	Left	Junior: 12.5-2	10.4" (27 cm)
28U33=R30-33	Right		
28U33=L33-36	Left	Junior: 2-4/4.5	11.6" (29.5 cm)
28U33=R33-36	Right		

### WalkOn Fit Kits

The *WalkOn Fit Kit* is an accessory that serves as test equipment to help you select the correct *WalkOn* ankle foot orthosis. It is used to assess whether the respective *WalkOn* orthosis is suitable for the patient with respect to the selected size and applicable indication. The following *Fit Kits* are available for the *WalkOn Lateral*:

Article number	Description	Four sizes included (adapted to the smallest size in each case)
28T6	28U33 WalkOn Lateral Fit Kit	R36-39, R42-45, L36-39, L42-45
28T16	28U33 WalkOn Lateral Junior Fit Kit	R24-27, R27-30, R30-33, R33-36 L24-27, L27-30, L30-33, L33-36