

28U33

WalkOn Lateral.

The *WalkOn Lateral* lifts the foot and also stabilizes the ankle joint.

Indication


Drop foot, with or without mild spasticity, e.g. caused by:

- Peroneal paralysis
- Stroke
- Traumatic brain injury
- Multiple sclerosis
- Neuromuscular atrophy

The indication must be determined by the doctor.

Benefits/advantages

- High energy return of the orthosis
- Effectively supports a physiological gait, even at different speeds and with different loads
- Stabilizes the ankle joint
- Increases the stride length by improving muscle dynamics
- The orthosis allows weight to be placed on the forefoot for climbing stairs and squatting
- Limits plantar flexion and supination when worn with a sturdy shoe
- Slim, lightweight design
- Optimized compliance thanks to a removable and washable pad/strap combination that can be aligned in the preferred direction for tightening the strap
- Optimized durability: adult sizes tested for two million cycles, or two million steps
- Material: carbon fiber prepreg
- Available in a variety of adult sizes





28U33 WalkOn Lateral

| Article number | Side | Shoe size (US) | Calf height |
|----------------|-------|---------------------------|-----------------|
| 28U33=L36-39 | Left | Men: 5-6 Women: 6-8.5 | 12.8" (32.5 cm) |
| 28U33=R36-39 | Right | | |
| 28U33=L39-42 | Left | Men: 6-9 Women: 8.5-10 | 13.8" (35 cm) |
| 28U33=R39-42 | Right | | |
| 28U33=L42-45 | Left | Men: 9-12 Women: 10-13 | 14.6" (37 cm) |
| 28U33=R42-45 | Right | | |
| 28U33=L45-48 | Left | Men: 12-15 | 15.4" (39 cm) |
| 28U33=R45-48 | Right | | |

WalkOn Fit Kits

The *WalkOn Fit Kit* is an accessory that serves as test equipment to help you select the correct *WalkOn* ankle foot orthosis. It is used to assess whether the respective *WalkOn* orthosis is suitable for the patient with respect to the selected size and applicable indication. The following Fit Kit is available for the *WalkOn Lateral*:

| Article number | Description | Four sizes included (adapted to the smallest size in each case) |
|----------------|------------------------------|--|
| 28T6 | 28U33 WalkOn Lateral Fit Kit | R36-39, R42-45, L36-39, L42-45 |