

WillowWood One Accessories

LIMBLOGIC® EXPULSION VALVE



Composed of a base, a one-way valve, and a release plug, the LimbLogic Expulsion Valve creates an easier donning and doffing process for an amputee using elevated vacuum. The valve, designed for transfemoral applications, is to be used in laminated sockets.

PRODUCT WEIGHT	15 g
PATIENT WEIGHT LIMIT*	350 lb (160 kg) for U.S. Activity Levels 1-3 300 lb (135 kg) for U.S. Activity Level 4

WARRANTY	12 months from date of patient fitting
-----------------	--

ITEM NUMBERS

Product	Item #
LimbLogic® Expulsion Valve	LLV-EV

ALIGNMENT BUTTONS **



Alignment buttons aid in consistent donning and orientation of the WillowWood One TF System. The buttons provide a physical lock between the removable brim and the definitive socket.

PATIENT WEIGHT LIMIT*	350 lb (160 kg) for U.S. Activity Level 3; 300 lb (136 kg) for U.S. Activity Level 4
------------------------------	---

ITEM NUMBERS

Product	Item #
Alignment Buttons, Replacement, Qty: 3	ONE-AB3

ONE SEAL **



The One Seal features 'fins' that create an airtight seal between the liner and the definitive socket. As part of the WillowWood One TF system, the silicone One Seal accommodates limb motion while maintaining a secure, airtight socket.

WARRANTY	12 months from date of patient fitting
-----------------	--

ITEM NUMBERS

Product	Item #
One Seal, Replacement, Medium Plus	ONE-SEAL-MP
One Seal, Replacement, Large	ONE-SEAL-LG
One Seal, Replacement, Large Plus	ONE-SEAL-LP
One Seal, Replacement, Extra Large	ONE-SEAL-XL

* Body weight plus any loads normally or routinely carried cannot exceed this weight limit.

** Requires WillowWood One certification to order.