

HOW TO GET CERTIFIED

1 TRAINING COURSE AT WILLOWWOOD'S EDUCATION CENTER

This one-day seminar covers how to fit patients with the WillowWood One® System. Course provides hands-on experience in critical aspects of the system process, such as evaluating the diagnostic socket fittings and customizing the sealing components.

2 LIVE EDUCATIONAL WEBINAR

This online, instructor-led training session (via WebEx) provides an overview of key aspects of fitting patients with the WillowWood One System. Ample opportunity for interaction with the clinician instructor is available.

FOR MORE INFORMATION ON HOW TO GET CERTIFIED, CALL **877.665.5443** OR VISIT **WILLOWWOODCO.COM**

HOW TO ORDER

1 Get certified using one of WillowWood's approved methods.

2 Order the WillowWood One TF System from WillowWood or a distributor.
a. For a complete vacuum system, determine WillowWood One part number.

ONE-TFXX-XX-X

Model:	Suspension:	Liner Size:	Liner Type:
ONE = WillowWood One System	TFVD* = Transfemoral Vacuum Suspension - Distal TFVS* = Transfemoral Vacuum Suspension - Side	MP = Medium Plus LG = Large LP = Large Plus XL = Extra Large	S = Alpha SmartTemp® D = Alpha Duo®

* Vacuum part numbers include definitive socket fabrication.

b. For a complete suction system, determine liner size, then order corresponding liner and seal and OFAB-ONE-TFP

3 Create a diagnostic socket, fit, and assess for your patient.


4 Using the included prepaid shipping label, send diagnostic socket and completed Custom Solutions Order Form to WillowWood for definitive socket fabrication.

5 Receive definitive socket.

6 Deliver socket to patient and review use and care instructions.

SYSTEM COMPONENTS

- Definitive Socket
- Alpha® Liner/One Gel Sock/One Seal
- LimbLogic® Side Mount/Distal Mount* (optional)
- Removable Brim
- Alignment Buttons
- Expulsion Valve



WillowWood One TF System Separate/Replacement Part Numbers	
Description	Item Number
Fabricless Alpha SmartTemp with One Gel Sock; Medium Plus	T990-1944
Fabricless Alpha SmartTemp with One Gel Sock; Large	T990-1946
Fabricless Alpha SmartTemp with One Gel Sock; Large Plus	T990-1948
Fabricless Alpha SmartTemp with One Gel Sock; Extra Large	T990-1947
Alpha Duo Liner with One Gel Sock; Medium Plus	D990-1944
Alpha Duo Liner with One Gel Sock; Large	D990-1946
Alpha Duo Liner with One Gel Sock; Large Plus	D990-1948
Alpha Duo Liner with One Gel Sock; Extra Large	D990-1947
One Seal; Medium Plus	ONE-SEAL-MP
One Seal; Large	ONE-SEAL-LG
One Seal; Large Plus	ONE-SEAL-LP
One Seal; Extra Large	ONE-SEAL-XL
One Gel Sock; Medium Plus	ONE-SOCK-MP
One Gel Sock; Large	ONE-SOCK-LG
One Gel Sock; Large Plus	ONE-SOCK-LP
One Gel Sock; Extra Large	ONE-SOCK-XL
One System Definitive Socket Fabrication	OFAB-ONE-TFP
One System Alignment Replacement Buttons (Qty: 3)	ONE-AB3
LimbLogic Expulsion Valve	LLV-EV

LINER SIZE CHART FOR USE WITH WILLOWWOOD ONE

Alpha Liner size determines the size of the One Gel Sock and One Seal.

	Medium Plus	Large	Large Plus	Extra Large
Proximal Measurement	30 - 54 cm (12 - 21 in)	33 - 64 cm (13 - 25 in)	41 - 80 cm (16 - 31.5 in)	46 - 90 cm (18 - 35.5 in)
Distal Measurement	20 - 29 cm (8 - 11.5 in)	28 - 36 cm (11 - 14 in)	32 - 39 cm (12.5 - 15.5 in)	37 - 45 cm (14.5 - 18 in)

WillowWood One TF U.S. L-Code Suggestions		
L Code	Product	Description
L5781* (vacuum only)	LimbLogic	Vacuum pump
L5652 (suction only)	Suction Suspension	Above knee suction suspension
L5679	Alpha Liner	Fabricless cushion liner
L8417	One Gel Sock	Gel sheath/sock
L5701	One Socket	Above knee socket
L5651	Removable Brim	Flexible inner
L5950	One Socket	Ultra-light material
L5920	One Socket	Alignable system
L5649	One Socket	Ischial containment
L5650	One Socket	Total contact
L5631	One Socket	Acrylic socket

* Only applies to WillowWood One system using elevated vacuum.

HOW TO ORDER

Certification is not required to purchase the WillowWood One TT System. Attending an education course on the WillowWood One System is highly recommended. All WillowWood One education offerings cover both transfemoral and transtibial applications.

- 1
- Order the WillowWood One TT System from WillowWood or a distributor.
a. Refer to willowwoodco.com to verify your patient is a candidate for the system.
b. Determine WillowWood One System part number.

ONE-TTXX-XX-XX


Model:	Suspension:	Liner Size:	Sleeve Size:
ONE = WillowWood One System	TTVD* = Transtibial Vacuum Suspension TTSU† = Transtibial Suction Suspension	MD = Medium MP = Medium Plus LG = Large LP = Large Plus XL = Extra Large	SM = Small MD = Medium LG = Large XL = Extra Large

* Definitive socket fabrication included.
† To include definitive suction socket fabrication, order OFAB-ONE-TTP

- 2
- Create a diagnostic socket, fit, and assess your patient.
- 3
- Two options for fabrication of the definitive socket:
a. Using the pre-paid shipping label, send diagnostic socket and completed order form to WillowWood for definitive socket fabrication. WillowWood will ship the definitive socket upon completion.
b. Fabricate the definitive socket at your facility. Please refer to and follow the provided fabrications instructions or watch the fabrication video at willowwoodco.com.
- 4
- Deliver system to patient and review use and care instructions.

ALPHA DUO® LINER SIZING CHART

Alpha Duo Liner size determines the size of the One Gel Sock.




SYSTEM KIT COMPONENTS

- Definitive Socket
- Alpha Duo® Liner/One Gel Sock/One Sleeve
- LimbLogic® Distal Mount* or Suction Pyramid (700-SP471)
- Flexible Inner Socket
- One Link



ONE SLEEVE SIZING CHART



WillowWood One TT Replacement Part Numbers

Description	Item Number
Alpha Duo Liner with One Gel Sock; Medium	D990-1940
Alpha Duo Liner with One Gel Sock; Medium Plus	D990-1944
Alpha Duo Liner with One Gel Sock; Large	D990-1946
Alpha Duo Liner with One Gel Sock; Large Plus	D990-1948
Alpha Duo Liner with One Gel Sock; Extra Large	D990-1947
One Gel Sock; Medium	ONE-SOCK-MD
One Gel Sock; Medium Plus	ONE-SOCK-MP
One Gel Sock; Large	ONE-SOCK-LG
One Gel Sock; Large Plus	ONE-SOCK-LP
One Gel Sock; Extra Large	ONE-SOCK-XL
One Sleeve; Small	ONE-SLEEVE-TT-SM
One Sleeve; Medium	ONE-SLEEVE-TT-MD
One Sleeve; Large	ONE-SLEEVE-TT-LG
One Sleeve; Extra Large	ONE-SLEEVE-TT-XL
One Link Kit	ONE-TT-LINKKIT

WillowWood One TT Fabrication Options

Description	Item Number
Fabrication Kit w/One Link & Sleeve Dummy; Small	ONE-TT-FABKIT-SM**
Fabrication Kit w/One Link & Sleeve Dummy; Medium	ONE-TT-FABKIT-MD**
Fabrication Kit w/One Link & Sleeve Dummy; Large	ONE-TT-FABKIT-LG**
Fabrication Kit w/One Link & Sleeve Dummy; Extra Large	ONE-TT-FABKIT-XL**
One System Transtibial Definitive Socket Fabrication	OFAB-ONE-TTP

** Sleeve size determines final digits in part number.

WillowWood One TT U.S. L-Code Suggestions

L Code	Product	Description
L5781* (vacuum only)	LimbLogic	Vacuum pump
L5647 (suction only)	Suction Suspension	Below knee suction suspension
L5679	Alpha Liner	Fabricless cushion liner
L8417	One Gel Sock	Gel sheath/sock
L5700	One Socket	Below knee socket
L5645	Flexible Inner Socket	TT Flexible inner socket
L5940	One Socket	Ultra-light material
L5910	One Socket	Alignable system
L5637	One Socket	Total contact
L5629	One Socket	Acrylic socket
L5685	One Sleeve	Suspension/sealing sleeve

* Only applies to WillowWood One system using elevated vacuum.