



CINDY WRIST SUPPORT

Application Instructions

Hand Wash Cold • Air Dry



1. Open brace. Separate hook closure. Place hand and forearm on brace.



2. Attach thumb strap to brace.



3. Slip hook and loop closure strap through ring and fasten securely.



4. Adjust for comfort.



Single Patient Use

BIRD & CRONIN INC.

Quality Medical Products To Help Improve Your Life

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail sales@birdcronin.com • www.birdcronin.com

©2012 Bird & Cronin, Inc. All Rights Reserved - 10005541A