

circaid reduction kit whole leg

circaid reduction kit - lower leg										
		Circumference in cm	regular				wide			
Step 1: Select lower leg width		cC	< 65				< 80			
		cB	< 55				< 65			
Step 2: Select lower leg length		Length in cm	short	standard	long		standard	long		
		ℓ mall - KC	< 32	32 - 36	> 36		≤ 36	> 36		
			▼				▼			
circaid reduction kit - upper leg										
		Circumference in cm	regular		wide			wide		
Step 3: Select upper leg width		cG	< 80		< 125			< 125		
		cE2	< 65		< 110			< 110		
Step 4: Select upper leg length		Length in cm	short	standard	short	standard	long	short	standard	long
		ℓ KC - Groin	< 37	≥ 37	< 37	37 - 41	> 41	< 37	37 - 41	> 41

