


Sizing

Inside Shoe Dimensions (in Inches)

Choice of shoe size is significantly influenced by foot swelling and volume of bandages/cast worn in conjunction with a DARCO shoe. To ensure the best fit, measure the widths and length your injured foot (including any dressings) and select the best shoe size using the graphic below.



US Shoe Size Reference Table

Use the following sizing table as a reference only.

Men's

Size	Shoe size	Catalog Number
Small	6-8	MSM1N
Medium	8.5 – 10	MSM2N
Large	10.5 – 12	MSM3N
X-Large	12.5 – 14	MSM4N

Women's

Size	Shoe size	Catalog Number
Small	4 – 6	MSW1N
Medium	6.5 – 8	MSW2N
Large	8.5 – 10	MSW3N

Children's

Size	Shoe size	Catalog Number
Pediatric	12-1	MSPN