

Fillauer Reciprocating Gait Orthosis (RGO)

Practitioner _____ Branch Name _____
PO# _____ Date: _____ Ship Via: _____
Address _____
City _____ State _____ Zip _____



Patient Information

Patient's Name _____
Age _____ Sex _____ Weight _____ Height _____
Ship To Address _____

☐ Fabricate Complete RGO ☐ Fabricate Pelvic Section Only

Measurements

1. Provide all measurements, especially the M-L diameters taken with caliper gauge

2. Supply a profile tracing from axilla to feet. Please check that tracing matches calipered diameters.

3. Plaster molds of each extremity are required.
A. Maintain ankle joint in neutral position for low heel shoe.
B. Correct ankle inversion or eversion as much as possible.
C. Reduce knee valgus or varus.
D. Use a right angle casting board.

4. When a plastic pelvic band is desired, refer to special casting suggestions in the RGO Fabrication Manual.

M-L at Metatarsal Head

RGO Options

Cable Type
☐ Horizontal Cable ☐ Dual Cable
☐ Rocker Bar (Metal Pelvic Section Only)

Pelvic Section
☐ Plastic Band ☐ Metal Band

Uprights
☐ Single ☐ Double

Hip Joints
☐ Regular Thrust Bearing Hip Joints
☐ Abduction Hip Joints
☐ Non-Abduction Hip Joints

Hip Joint Sizes
☐ Small (3/16x5/8)
☐ Medium (1/4x13/16)
☐ Large (5/16x7/8)

Knee Joint
☐ Fillauer Cam Lock
☐ Drop Locks ☐ Bail Lock
☐ Other _____

Knee Joint Bar Size
☐ 3/16 x 5/8
☐ 3/16 x 3/4 ☐ 1/4 x 3/4

Heel Height
☐ 0" ☐ 1/4" ☐ 1/2"

Additional Instructions

