#### compression levels available: 20-30, 30-40 mmhg

mediven® harmony sleeves, gauntlets and gloves - 20-30 and 30-40 mmhg size chart

circumference	1	2	3	4	5	6	7	8
"g" extra-wide	28-31	30-34	32.5-36.5	34.5-39.5	36.5-41.5	39-44	41-46.5	44.5-50.5
"g" standard	25-28	27.5-31	29.5-33.5	31.5-36	33-37.5	35.5-40.5	37.5-42.5	40-45.5
"e" extra-wide	24-26.5	26-29.5	28-31.5	29.5-34	31-35.5	33-37.5	35-40	38-43
"e" standard	23-25	24-27	25.5-29	27-31	28.5-32.5	30-34.5	32-36.5	34.5-39.5
"c1"	15-16.5	16-18.5	17.5-20	19-22	20.5-24	22.5-26	24.5-28.5	27-31
"c"	13.5-14.5	14-16	15.5-18	17-20	18.5-21.5	20-23.5	22-25.5	24.5-28.5
"a"	15-18	15.5-19.5	16.5-20.5	18.5-23	19.5-24	21.5-26	23.5-28	25.5-30

gloves & gauntlets available in sizes 2-7

arms leeves available in sizes 1-8

#### **MEASURING INSTRUCTIONS**

- 1. Select product/style from measuring chart.
- 2. Measure at indicated locations on diagram.
- 3. Locate proper product size based on measurements.

## For All Arm Sleeves:


Measure **G e g** 


### For All Gauntlets & Gloves:

Measure (a) and (C)

# For All One-Piece Arm Sleeves & Gauntlets:

Measure a ceg





Circumference measurements for gloves and gauntlets.