

FLEX-FOOT® JUNIOR



Flex-Foot Junior provides that important energy return that children require. The Carbon X® Active Heel stores energy and absorbs shock loads while the full-length toe lever contributes to stability and even stride length.

The low profile, energy storing foot module features a sandal toe shape and an integrated male pyramid. The shape of the pyramid adapter contributes to progressive stiffening of the foot, and the carbon fiber layering technique is specifically designed to offer the support and flexibility needed for the varied movements of children.



USER INFORMATION	
Amputation Level:	Transfemoral and Transtibial
Impact Level:	Low to High
Maximum Patient Weight:	45kg (100lbs)

FOOT INFORMATION	
Categories:	1-3
Sizes:	16-21
Weight of Foot: (Size 19)	288g (10.2oz)
Build Height: (Size 19)	46mm (1 13/16")
Heel Height:	10mm (3/8")
Adapter Options:	Integrated Male Pyramid




CATEGORY SELECTION CHART

WEIGHT KG	15-24	25-34	35-45
WEIGHT LBS	33-53	54-75	76-99
CATEGORY	1	2	3

NOTE: The above categories apply only to Flex-Foot Junior
*All sizes not available for each category



CLEARANCE



FLEX-FOOT® JUNIOR

KITS

FLEX-FOOT JUNIOR KIT WITH MALE PYRAMID




Includes:

- Toe Cushion
- Foot Module
- Heel Module
- Heel Fasteners
- Heel Wedges
- Male Pyramid

COMPONENTS

JUNIOR HEEL WEDGE KIT


PART#	FOOT SIZE	CAT.
CM150017	16-21	ALL



Includes:

- 3 Wedges
- Adhesive

COSMETIC FINISHING



Includes:

- Foot Cover*
- Flex Foot Sock

*Please specify when ordering brown covers add "BR" to the part#.

PART#	DESCRIPTION	SIZE
FCX63008	Flex-Foot Sock	16-21

