

Philadelphia Adjustable Collar

Part#	Size/Description
PHP-AA	Conforms to sizes 2, 3 and 4 of the original Philadelphia Tracheotomy Collar

NOTE: This product is recommended for single patient use. Do not reuse on multiple patients. Follow a qualified healthcare provider's instructions regarding usage and infection control.

Philadelphia Adjustable Collar - Sizing reference

Neck height	
The neck height adjustment accommodates 2 ¼" (6 cm), 3 ¼" (8 cm) and 4 ¼" (11 cm)	
Circumference size	Range
Small (2 tears)	10"-13" (25.4-33.0 cm)
Medium (1 tear)	13"-16" (33.0-40.6 cm)
Large (no tear)	16"-19" (40.6-48.3 cm)

NOTE: The table above represents common sizing only. Sizing must be verified by the healthcare provider for each patient individually.