

The Boxer™ allows the caregiver an option to casting. Meticulous attention to every aspect of the Boxer™ is apparent. The Palmar Spoon, Ulnar/Dorsal, and Metacarpal Stays are designed to eliminate pressure on the ulnar styloid and hold the wrist and fingers in the preferred angle of support when treating 4th and 5th metacarpal fractures.

SOXER

U.S. & FOREIGN PATENTS PENDING.

Materials:

Orthopaedic Neoprene

Fit & Function:

- Comfortable
- Alternative to casting or splinting
- Can be applied with one hand
- Allows adjustment for swelling and removal for skin care
- Velcro® Hoop Closure will not snag clothing
- Removable Palmar Spoon, Ulnar/Dorsal, and Metacarpal Stays
- Ulnar/Dorsal Stay is pre-adjusted with fingers in 70° flexion
- Palmar Stay is pre-adjusted with wrist in 30° extension
- Metacarpal Stay is easily adjustable
- Engineered Dorsal Stay eliminates ulnar styloid pressure
- Semi-Universal

Options & Ordering: 11357 Boxer™, Left 11358 Boxer™, Right Measure the circumference of the affected forearm 7 inches above the wrist joint. Pediatric/Youth size

determined by body weight.		
Pediatric/Youth	30 to 50 Lbs.	13.5 to 22.5 kg
XSmall/Small	8 to 10 inches	20 to 25 cm
Medium/Large	10 to 12 inches	25 to 30.5 cm
XLarge/XXLarge	11.5 to 15 inches	29 to 38 cm

Suggested HCPCS Codes:

- L3917 Hand orthosis, metacarpal fracture orthosis, prefabricated, includes fitting and adjustment.
- L3807 Wrist hand finger orthosis, without joint(s), prefabricated, includes fitting and adjustments, any type.



PROFESSIONAL PRODUCTS, INC. 54 HUGH ADAMS ROAD DEFUNIAK SPRINGS, FL 32435

USA & CANADA (800) 234-9004 INTERNATIONAL (850) 892-5731 For more information on **EZYWRAP**® products and services, contact your