




**EXotec®**  
**CUBITAL  
TUNNEL  
SPLINT**



Cubital tunnel syndrome  
is the second leading  
compression neuropathy...





Cubital tunnel syndrome occurs when the ulnar nerve becomes trapped or irritated as it passes through the cubital tunnel. For functional and comfortable splinting, our Cubital Tunnel Splint positions the elbow to minimize movement and reduce irritation (especially at night).

**Materials:**

- Orthopaedic neoprene with Kydex® ventral insert

**Fit & Function:**

- Superior design for fit and comfort
- Easy application
- Adjustable for variable pre or post-op protocols
- Removable option to casting
- Allows for skin care and edema relief
- Fits left and right
- Universal size

**Ordering:**

11377-00 Cubital Tunnel Splint

**Suggested HCPCS Code:**

- L3762 Elbow orthosis, rigid, without joints, includes soft interface material, prefabricated, includes fitting and adjustment

