

### CARE:

To clean, unroll bandage from its core and rinse in luke warm water by hand. To prolong the life of the bandage's self-adhesive property, do not use detergent or scrub the bandage. Bandage will stick to itself through several wearings. After this it is necessary to use a fresh bandage. Washing will generally reduce the degree of self-adhesion.

### CAUTION:

This product contains natural rubber latex which may cause allergic reactions.

Fiber content: cotton, polyester, rubber

**CHAMPION®**  
**championsupports.com**

### 2" BANDAGE IS RECOMMENDED FOR:



**ANKLE**



**FOOT**



**WRIST**

- Inconspicuous neutral color.
- Can be washed and reused.
- Clips allow bandage to be firmly secured.

### CARE:

To clean, unroll bandage from its core and rinse in lukewarm water by hand. To prolong the life of the bandage's self-adhesive property, do not use detergent or scrub the bandage. Bandage will stick to itself through several wearings. After this it is necessary to use a fresh bandage. Washing will generally reduce the degree of self-adhesion.

### CAUTION:

This product contains natural rubber latex which may cause allergic reactions.

Fiber content: cotton, polyester, rubber

### 4" BANDAGE IS RECOMMENDED FOR:



**KNEE**



**SHOULDER**



**LOWER LEG**



**UPPER LEG**

- Inconspicuous neutral color.
- Can be washed and reused.
- Clips allow bandage to be firmly secured.

**CHAMPION®**

**championsupports.com**

### CARE:

To clean, unroll bandage from its core and rinse in lukewarm water by hand. To prolong the life of the bandage's self-adhesive property, do not use detergent or scrub the bandage. Bandage will stick to itself through several wearings. After this it is necessary to use a fresh bandage. Washing will generally reduce the degree of self-adhesion.

### CAUTION:

This product contains natural rubber latex which may cause allergic reactions.

Fiber content: cotton, polyester, rubber

### 3" BANDAGE IS RECOMMENDED FOR:



**ANKLE**



**FOOT**



**LOWER LEG**



**UPPER LEG**

- Inconspicuous neutral color.
- Can be washed and reused.
- Clips allow bandage to be firmly secured.

**CHAMPION®**  
**championsupports.com**