

A man in a dark long-sleeved shirt, tan cargo shorts, and a black cap is in the middle of a golf swing. He is wearing a black knee brace on his left knee. The background shows autumn foliage and a clear sky.

UNLOAD
THE PAIN

GET YOUR PATIENTS
BACK IN THE SWING

OA NÖVEL KNEE BRACE



Experience the United Ortho Difference.

UNITED
ORTHO

OA NOVEL

FEATURING THE DYNAMIC UNLOADING HINGE

How It Works: The Dynamic Unloading Hinge lengthens the strut on the affected side, (in the last 30 degrees of extension), which changes the angle of the brace and unloads the knee—and the pain!





OA NōVel Knee Braces



Cinching Dial
makes incremental
adjustments a cinch.

**Dynamic Unloading
Hinge** unloads the
knee—and the pain!

**Load Dispersion
Cuffs** provide
increased comfort.



Adjustment Dial
delivers increased
unloading, if needed.



**Quick Release
Buckles** for easy
on and off.

Experience the United Ortho Difference.

OA NōVEL MEDIAL AND LATERAL

Part #	Description
301261-XX	OA, NōVel, OTS, Medial, Right
301262-XX	OA, NōVel, OTS, Medial, Left
301271-XX	OA, NōVel, OTS, Lateral, Right
301272-XX	OA, NōVel, OTS, Lateral, Left

For accessories, visit www.united-ortho.com

Examples of Use

- Unicompartmental OA-medial and lateral
- Chondral defect procedures
- Meniscal procedures

One year warranty on frames and hinges.

HCPCS- L1845/K0902



OA NōVel - Medial



OA NōVel - Lateral



SIZING CHART

XX =	Size	Circ. Above*	Circ. Below**
01	XS	13.00"-15.50"	10.00"-11.75"
03	SM	15.50"-18.00"	11.75"-13.50"
05	MD	18.00"-20.50"	13.50"-15.50"
07	LG	20.50"-23.00"	15.50"-17.50"
08	XL	23.00"-26.00"	17.50"-19.25"
09	2XL	26.00"-29.00"	19.25"-21.50"

*Circumference measurement taken 6" above mid-patella.

**Circumference measurement taken 6" below mid-patella.



1-800-227-8748 | UNITED-ORTHO.COM

UNITED ORTHO ENTERPRISES

2235 PENNSYLVANIA STREET | FORT WAYNE, IN 46803 | P 800-227-8748 | F 877-688-7975